

NEW

MediKit

First Aid & Safety Products

BS8599-1:2019

Compliant First Aid Kits

WE'VE GOT
FIRST AID
COVERED

t.+44 (0)1522 513 512

e.sales@medikit.co.uk

www.medikit.co.uk

NEW

New BS8599-1:2019 Compliant First Aid Kits

MediKit

First Aid & Safety Products

Introducing the NEW BS8599-1:2019 Compliant First Aid Kits. The BSI workplace first aid kit (BS8599-1) content has been under review to ensure the contents are still suitable and fit for purpose. Modern risk and injuries have been taken into consideration by the BSI and other stakeholders such as the BHTA.

Following this review there have been some amendments to the existing kit contents and also the introduction of two new kits; a Personal Injury Kit and a Critical Injury Pack.

What's changed?

- Two new kits, a Personal Issue Kit and a Critical Injury Pack. Both kits will be BS8599-1:2019 compliant
- BS8599-1 Travel Kit has been replaced by BS8599-2 Medium Travel and Motoring Kit
- Minor changes to the workplace first aid kits Small, Medium and Large



What's new?

The Critical Injury Packs and the Personal Issue Kits are the first of their kind to meet BSI regulations; born from the concern of heightened risk for mass casualty situations. The Critical Injury Packs are designed for anywhere a person may suffer with a severe injury. The kits would be ideal for people working in manufacturing and machinery or construction and forestry. The content included will help to manage severe blood loss. We recommend conducting a risk assessment to decide the appropriate amounts of Critical Injury Packs are provided for both employees and the public in that location at any one time.

When do the New BS8599-1:2019 Compliant First Aid Kits take effect?

Manufacturers have a year to update their products to the new standard meaning any existing stock will still meet the current British Standard. Medikit have already implemented these changes and switched over to the new standard to best support our distributors.

Questions? Get in touch our team are always happy to help

With over 25 years experience we operate under the ISO:9001 Quality Accreditation and all our first aid products carry relevant CE certification. Contact our team today for any questions.



Standard Range

Compliant with BS8599-1:2019

- Tough and durable plastic box
- Optional clip in wall bracket

32406	Standard First Aid Kit - Small Content Plastic Box - 240 x 170 x 80mm
32408	Standard First Aid Kit - Medium Content Plastic Box - 270 x 220 x 90mm
32410	Standard First Aid Kit - Large Content Plastic Box - 340 x 270 x 100mm



Premium Range

Compliant with BS8599-1:2019

- Stylish plastic box complete with internal compartments
- Complete with bracket

34406	Premium First Aid Kit - Small Content Plastic Box & Wall Bracket - 280 x 280mm
34408	Premium First Aid Kit - Medium Content Plastic Box & Wall Bracket - 280 x 280mm
34410	Premium First Aid Kit - Large Content Plastic Box & Wall Bracket - 340 x 330mm



Deluxe Range

Compliant with BS8599-1:2019

- Robust ABS plastic boxes with internal compartments
- Rubber lid seal reducing contamination from dust and moisture
- Complete with clip in wall bracket

33406	Deluxe First Aid Kit - Small Content ABS Box & Wall Bracket - 270 x 185 x 115mm
33408	Deluxe First Aid Kit - Medium Content ABS Box & Wall Bracket - 320 x 220 x 130mm
33410	Deluxe First Aid Kit - Large Content ABS Box & Wall Bracket - 320 x 220 x 130mm



Guidance to choosing the right kit

Category of Hazard	Number of Employees	Size of First Aid Kit
Low Hazard Environment e.g Offices, Shops, Libraries etc	Less than 25	1 x Small Kit
	25 to 100	1 x Medium Kit
	More than 100	1 x Large Kit Per 100 Employees
High Hazard Environment e.g Engineering, Assembly Work, Food Processing, Construction etc	Less than 5	1 x Small Kit
	5 to 25	1 x Medium Kit
	More than 25	1 x Large Kit Per 25 Employees

Kit Content

BS8599-1:2019 Compliant	Small	Medium	Large
First Aid Guidance Leaflet	1	1	1
Medium Sterile Dressing 12 x 12cm	2	4	6
Large Sterile Dressing 18 x 18cm	2	3	4
No.16 Sterile Eye Pad	2	3	4
Triangular Bandage	2	3	4
Adhesive Plasters Assorted Sizes*	40	60	100
Moist Cleansing Wipes	20	30	40
Disposable Nitrile Glove	6	9	12
Burn Dressing 10 x 10cm	1	2	2
Sterile Finger Bandage	2	3	4
Microporous Tape	1	2	3
Resuscitator Face Shield	1	1	2
Emergency Foil Blanket	1	2	3
Heavy Duty Clothing Shears	1	1	1
Conforming Bandage 7.5cm x 4m	1	2	2
Sterile Eye Wash Solution 250ml	-	-	-

*Blue Detectable Plasters in kits for Food Processing & Catering Areas

Travel and Motoring Kit

Compliant with BS8599-1:2019

- Choice of boxes with optional wall brackets or grab bag for easy stowage within a vehicle

32464	Travel and Motoring Kit Box c/w Bracket - 240 x 170 x 80mm
32465	Travel and Motoring Kit Pouch - 210 x 140 x 60mm
32466	Travel and Motoring Kit Deluxe Pouch - 255 x 155 x 95mm



BS8599-1:2019 & BS8599-2:2014 Compliant

HSE Guidance Leaflet	1
Contents Label	1
Sterile 12 x 12 cm Medium Dressing	1
Washproof Plasters (Pack 10)	1
Triangular Bandage	1
Nitrile Gloves (Pair)	2
Moist Cleansing Wipes	10
Large Adhesive Dressing 7cm x 6cm	1
Trauma Dressing 10 x 18cm	1
Burn Dressing 10 x 10cm	2
Foil Blanket	1
Heavy Duty Clothing Shears	1
Resusciade	1

Personal Issue Kit

Compliant with BS8599-1:2019

- Choice of boxes and bags available to suit your needs assessment

31110	Personal Issue Kit Box - 227 x 165 x 47mm
31112	Personal Issue Kit Pouch - 180 x 110 x 45mm
31114	Personal Issue Kit Bum Bag - 420 x 125 x 120mm



BS8599-1:2019 Compliant

Guidance Leaflet	1
Contents Label	1
Large 18cm x 18cm Dressing	1
Triangular Bandage	1
Washproof Plasters Pack 10	1
Moist Cleansing Wipes	4
Nitrile Gloves (Pair)	2
Resusciade	1
Heavy Duty Clothing Shears	1
Foil Blanket	1

Critical Injury Pack

Compliant with BS8599-1:2019

- Available in either Standard or Premium
- Designed for 'drop or throw' in areas with mass casualties

31120	Critical Injury Pack Standard - 230 x 150 x 85mm
31122	Critical Injury Pack Premium - 230 x 150 x 85mm

Our BS8599-1:2019 compliant critical injury pack provides supplies for treatment of severe bleeds and trauma. This kit is the first of its kind to be compliant with the latest BSI standards; born from the concern of heightened risk for mass casualty situations.

Supplied complete with content, the kits include products to help manage severe blood loss such as a tourniquet and haemostatic dressing.

Choose from our Standard or Premium Pack. Premium packs include Celox™ Gauze Z-Fold haemostatic dressing because they have been tested in controlled environments and in action. The results show you can save lives with the use of Celox™. The dressings have Celox™ granules bonded onto the gauze, therefore giving the dressing the same effective technology as the original Celox™ Gauze roll.

The Premium pack also includes a military CAT tourniquet. The US Army use these tourniquets and deem it to be 100% effective. Since its introduction in 2005, the CAT tourniquet has dropped mortality rate due to extreme blood loss by 85%.* The tourniquet is a true one handed device, simply apply to the affected area, wind to stop blood flow and secure in place with the velcro strap.

BS8599-1:2019 Compliant Standard Premium

Guidance Leaflet	1	1
Contents Label	1	1
Nitrile Gloves (Pair)	2	2
Heavy Duty Clothing Shears	1	1
Large Trauma Dressing 15cm x 18cm	2	2
WoundClot Haemostatic Wound Dressing	2	-
Celox Haemostatic Wound Dressing	-	2
Tourniquet	1	-
CAT Military Tourniquet	-	1
Foil Blanket	1	1

Our standard pack is suitable for areas deemed as low risk. The kit includes standard haemostatic dressings and also a quick release tourniquet. Our quick release tourniquet is a buckle system, simply apply around the area and pull tight to stem blood flow. Easy to remove, simply press the button on the buckle for quick release.

We recommend carrying out an assessment of needs because this will ensure you provide the appropriate amount of Critical Injury Packs. Both employees and the public should be a consideration for that location at any one time.



*figures quoted at the time of publishing by CAT Resources LLC