

MediKit

First Aid & Safety Products



MediKit Ltd has provided a quality range of First Aid Products with excellent service to our customers since 1996. Our customers benefit from next day delivery on orders for our core range, which includes:

First Aid Kits, Eye Wash Stations, Plasters and Consumable Dressings, plus an expanding range of First Aid Room and Evacuation and Rescue items.

MediKit is Quality Assured to ISO 9001 and our First Aid Dressings carry CE 0120 accreditation, providing total quality assurance on many of our products and a service our customers can depend on.

To Place An Order

Phone: +44 (0) 1522 513 512
Email: sales@medikit.co.uk
Web: www.medikit.co.uk

Payment Methods

BACS: Please contact our Sales Office for full bank details
Credit/Debit Card: Company or Personal cards
Cheque: Made payable to MediKit Limited

Terms and Conditions Summary

Payment: 30 days net
Title: Title in goods does not pass until payment is received in full
VAT: Currently charged at 20%
Carriage: Paid on orders over £75 nett order value, UK Mainland

Resale
Customer: Paid on orders over £150 nett order value
Small Order
Charge: £3.00 on orders under £30 nett order value
Minimum order value: £10 nett order value
Samples: Will be invoiced in all cases. A full credit will be issued, not including carriage costs, upon return of goods in a saleable condition. Sterile and limited shelf life products cannot be returned for credit without prior written authority from MediKit Limited.



www.medikit.co.uk

Visit our website for up to date company and product information. Digital versions of this catalogue are available online via the Downloads section.

Accreditation

MediKit Limited's high quality range of first aid sterile bandages and plasters are accredited by the notified body SGS Limited and carry the number 0120.

Risk Assessment

The Health & Safety (First Aid) Regulation 1981 and the 1997 revised Approved Code of Practice.

The above states an employer shall provide, or ensure that there are provided, such equipment and facilities which are adequate and appropriate in the circumstances for enabling first aid to be rendered to employees if they are injured or become ill at work.

Please Note: All specifications throughout this catalogue are approximate and may be subject to change without prior notification.

MediKit© All material copyright and not to be reprinted, reproduced or stored in a retrieval system without prior written permission obtainable from MediKit Limited, Lincoln, LN2 4SY.

**VISIT US
ONLINE**

Workplace First Aid Kits [Compliant with BS8599-1]	4-5
Motor Vehicle First Aid Kits [Compliant with BS8599-2]	6
First Aid Kits [Compliant with BS8599-1]	7
HSE Travel and Lone Worker Kits	8
HSE Workplace First Aid Kits	9
Empty Metal Cabinets	10
Metal Cabinet First Aid Kits	11
Travel and Domestic First Aid Kits	12
Commercial Vehicle First Aid Kits	13
Incident and Evacuation First Aid Kits	14-17
Forestry First Aid Kits	18
Farming and First Response First Aid Kits	19
Burns First Aid	20-21
Sports First Aid Kits	22-23
Sports First Aid Consumables	24-25
Plasters	26-27
Bandages	28-29
Dressings and Splints	30-31
Eyewash Stations and Solution	32-33
Bio Hazard Clean Up	34
Sharps Handling and Disposal	35
First Aid Room Furniture Packages	36-37
First Aid Room	38-39
Evacuation and Rescue Equipment	40-43
Defibrillators and Accessories	44-45
First Aid Equipment and Accessories	46-47
Workplace Fire Extinguishers	48-49
Vehicle Fire Extinguishers and Safety	50-51
Fire Safety Products	52-53
First Aid Boxes and Bags	54-57
First Aid Training Manikins	58
Cotec Finger Protection Tape	59
Disposable First Aid Gloves	60
Signs	61

A-Z product listing can be found at the back of this catalogue.

Workplace First Aid Kits

Compliant with BS8599-1:2019

MediKit's range of BS8599-1:2019 Compliant Workplace First Aid Kits ensure employers are fully up to date with current regulations.

- Updated content to suit the modern workplace; contains a more comprehensive inventory compared to the HSE 10, 20 and 50 Person First Aid Kits
- Available in three sizes Small, Medium and Large - size required dependent on the workplace risk, see help guide below.
- Providing BS8599-1: 2019 Compliant First Aid Kits in the workplace is seen as best practice



Guidance to choosing the right kit

Category of Hazard	Number of Employees	Size of First Aid Kit
Low Hazard Environment e.g Offices, Shops, Libraries etc	Less than 25	1 x Small Kit
	25 to 100	1 x Medium Kit
	More than 100	1 x Large Kit Per 100 Employees
High Hazard Environment e.g Engineering, Assembly Work, Food Processing, Construction etc	Less than 25	1 x Small Kit
	25 to 100	1 x Medium Kit
	More than 25	1 x Large Kit Per 25 Employees

Kit Content

BS8599-1:2019 Kit Content	Small	Medium	Large
Guidance Leaflet	1	1	1
Contents Label	1	1	1
No.16 Eye Dressing	2	3	4
Sterile 12 x 12 cm Medium Dressing	2	4	6
Sterile 18 x 18 cm Large Dressing	2	3	4
Washproof Plasters (Pack of 20)	2	3	5
Triangular Bandage	2	3	4
Nitrile Disposable Glove – Pair	6	9	12
Moist Cleansing Wipes	20	30	40
Microporous Tape 2.5cm x 10m	1	2	3
Conforming Bandage 7.5cm x 4m	1	2	2
Quick Fix Sterile Finger Dressing	2	3	4
Burn Dressing 10 x 10cm	1	2	2
Emergency Foil Blanket – 2.1m x 1.3m	1	2	3
Heavy Duty First Aid Snips	1	1	1
Resuscitator	1	1	2

*Blue Detectable Plasters in kits for Food Processing & Catering Areas.

Standard Range

Compliant with BS8599-1:2019

- Tough and durable plastic box
- Optional clip in wall bracket

32406	Standard First Aid Kit - Small Content Plastic Box - 240 x 170 x 80mm	A
32408	Standard First Aid Kit - Medium Content Plastic Box - 270 x 120 x 90mm	A
32410	Standard First Aid Kit - Large Content Plastic Box - 340 x 270 x 100mm	A



Premium Range

Compliant with BS8599-1:2019

- Stylish plastic box complete with internal compartments
- Complete with bracket

34406	Premium First Aid Kit - Small Content Plastic Box & Wall Bracket - 280 x 280mm	A
34408	Premium First Aid Kit - Medium Content Plastic Box & Wall Bracket - 280 x 280mm	A
34410	Premium First Aid Kit - Large Content Plastic Box & Wall Bracket - 340 x 330mm	A

Deluxe Range

Compliant with BS8599-1:2019

- Robust ABS plastic boxes with internal compartments
- Rubber lid seal reducing contamination from dust and moisture
- Complete with clip in wall bracket

33406	Deluxe First Aid Kit - Small Content ABS Box & Wall Bracket - 270 x 185 x 115mm	A
33408	Deluxe First Aid Kit - Medium Content ABS Box & Wall Bracket - 320 x 220 x 130mm	A
33410	Deluxe First Aid Kit - Large Content ABS Box & Wall Bracket - 320 x 220 x 130mm	A



Motoring First Aid Kits

Compliant with BS8599-2:2014

When employees are working away from the employer's establishment, the employer must still provide first aid cover e.g company cars and commercial vehicles.

BS8599-2 Compliant Motoring First Aid Kits are supplied in three sizes. Choice of kit is determined by the size of the vehicle and number of passengers.



Guidance to choosing the right kit

Type of Vehicle	Number of Passengers	Size of First Aid Kit
Moped, Motorcycle, Tricycle etc	Maximum of 3	Small
Car, Taxi, Van, Truck etc	Maximum of 8	Medium
Minibus, Bus, Coach	Maximum of 16	Large

Passengers exceeding 16 will require multiple kits

BS8599-2 Vehicle First Aid Kit Content	Small	Medium	Large
First Aid Guidance Leaflet	1	1	1
First Aid Kit Content List	1	1	1
Trauma Dressing Large 15 x 18cm	-	-	1
Trauma Dressing Medium 10 x 18cm	1	1	2
Medium Sterile Dressing 12 x 12cm	-	1	2
Triangular Bandage	-	1	2
Adhesive Dressing 7.5 x 6cm	-	1	2
Adhesive Plasters Assorted Sizes *	10	10	20
Moist Cleansing Wipes	5	10	20
Disposable Nitrile Glove	1	2	5
Burn Dressing 10 x 10cm	1	2	4
Resuscitator Face Shield	1	1	2
Emergency Foil Blanket	-	1	3
Heavy Duty Clothing Shears	1	1	1



Motor Vehicle First Aid Kits

Compliant with BS8599-2

32460	Motor Vehicle First Aid Kit - Small Content Plastic Box c/w Bracket 240 x 170 x 80mm	A
32462	Motor Vehicle First Aid Kit Small Content in Nylon Zip Pouch 180 x 110 x 45mm	A
32464	Motor Vehicle First Aid Kit - Medium Content Plastic Box c/w Bracket 240 x 170 x 80mm	A
32466	Motor Vehicle First Aid Kit Medium Content in Nylon Zip Pouch 255 x 155 x 95mm	A



32468	Motor Vehicle First Aid Kit - Large Content Plastic Box c/w Bracket 270 x 220 x 90mm	A
32470	Motor Vehicle First Aid Kit Large Content in Nylon Zip Pouch 255 x 155 x 95mm	A

Personal Issue Kit

Compliant with BS8599-1:2019

- Choice of boxes and bags available to suit your needs assessment

31110	Personal Issue Kit Box - 227 x 165 x 47mm	
31112	Personal Issue Kit Pouch - 180 x 110 x 45mm	
31114	Personal Issue Kit Bum Bag - 420 x 125 x 120mm	



BS8599-1:2019 Compliant

Guidance Leaflet	1
Contents Label	1
Large 18cm x 18cm Dressing	1
Triangular Bandage	1
Washproof Plasters Pack 10	1
Moist Cleansing Wipes	4
Nitrile Gloves (Pair)	2
Resusciate	1
Heavy Duty Clothing Shears	1
Foil Blanket	1

Critical Injury Pack

Compliant with BS8599-1:2019

- Available in either Standard or Premium
- Designed for 'drop or throw' in areas with mass casualties

31120	Critical Injury Pack Premium - 230 x 150 x 85mm	
31122	Critical Injury Pack Standard - 230 x 150 x 85mm	



Our BS8599-1:2019 Compliant critical injury pack provides supplies for treatment of severe bleeds and trauma. This kit is the first of its kind to be compliant with the latest BSI standards; born from concern of heightened risk of mass casualty situations.

Supplied complete with content, the kits include products to help manage severe blood loss such as tourniquet and haemostatic dressing.

BS8599-1:2019 Compliant Standard Premium

Guidance Leaflet	1	1
Contents Label	1	1
Nitrile Gloves (Pair)	2	2
Heavy Duty Clothing Shears	1	1
Large Trauma Dressing 15cm x 18cm	2	2
WoundClot Haemostatic Wound Dressing	2	-
Celox Haemostatic Wound Dressing	-	2
Tourniquet	1	-
CAT Military Tourniquet	-	1
Foil Blanket	1	1

HSE Travel & Lone Worker Kits

When employees are working away from the employer's establishment, the employer must still provide first aid cover e.g company cars and commercial vehicles.

HSE Standard Travel Kits

These kits are suitable for One Travelling Person and have the minimum HSE content.

Employers may wish to exceed this content and provide a BS8599-2 Compliant Motor Vehicle First Aid Kits. Please see page 6 for full details.



33111	HSE - Single Person First Aid Kit Bum Bag - 420 x 125 x 120mm	A
32114	HSE - Single Person First Aid Kit Plastic Box - 227 x 165 x 47mm	A
33115	HSE - Single Person First Aid Kit PVC Clear Zip Wallet - 260 x 210mm	A
33119	HSE - Single Person First Aid Kit PU Nylon Pouch with Belt Loop - 250 x 200mm	A

Contents	HSE 1 Person Travel	Lone Worker
First Aid Guidance Leaflet	1	1
Medium Sterile Dressing 12 x 12cm	-	2
Large Sterile Dressing 18 x 18cm	1	1
No.16 Sterile Eye Pad	-	1
Triangular Bandage	2	2
Adhesive Plasters Assorted Sizes	10	10
Moist Cleansing Wipes	4	6
Safety Pins	6	6
Disposable Nitrile Glove	1	1



Lone Worker Kits

These kits have the Standard HSE content for One Travelling Employee with additional items such as extra dressings and wipes. Available in plastic boxes or a belt mountable Nylon First Aid Pouch.

33103B	Lone Worker Kit - HSE Extra Content Plastic Box c/w Bracket Plastic Box - 210 x 140 x 52mm	A
32104D	Lone Worker Kit - HSE Extra Content Nylon Zip Pouch 180 x 110 x 45mm	A

HSE Workplace First Aid Kits

HSE First Aid Kits are being replaced in the workplace by BS8599-1 compliant kits, please see page 5 for more information. The HSE Standard First Aid Kits enable employers to comply affordably with the minimum requirements under current Health and Safety Regulations.

- Available in three sizes 10, 20 and 50 Person
- Cost effective refill kits available

Contents	1-10	1-20	1-50
First Aid Guidance Leaflet	1	1	1
Medium Sterile Dressing 12 x 12cm	6	9	12
Large Sterile Dressing 18 x 18cm	2	3	4
No.16 Sterile Eye Pad	2	4	6
Triangular Bandage	4	6	8
Adhesive Plasters Assorted Sizes*	20	40	60
Moist Cleansing Wipes	6	10	20
Safety Pins	6	12	12
Disposable Nitrile Glove	1	2	3

*Blue Detectable Plasters in kits for Food Processing



Deluxe Range

- Robust ABS plastic boxes with internal compartments
- Rubber lid seal reducing contamination from dust and moisture
- Complete with clip in wall bracket

B70107	HSE - Deluxe (1-10 person) ABS Box c/w Bracket - 270 x 185 x 115mm	A
B70108	HSE - Deluxe (1-20 person) ABS Box c/w Bracket - 320 x 220 x 130mm	A
B70112	HSE - Deluxe (1-50 person) ABS Box c/w Bracket - 320 x 220 x 130mm	A



Standard Range

- Tough and durable plastic box
- Optional clip in wall bracket

32106	HSE - Standard (1-10 person) Plastic Box - 240 x 170 x 80mm	A
32108	HSE - Standard (1-20 person) Plastic Box - 270 x 220 x 90mm	A
32110	HSE - Standard (1-50 person) Plastic Box - 340 x 270 x 100mm	A

When a wall bracket is required apply 'B' to the product code e.g 32106B



Premium Range

- Stylish plastic box complete with internal compartments
- Complete with bracket

34106	HSE - Premium (1-10 person) Plastic Box c/w Bracket - 280 x 280mm	A
34108	HSE - Premium (1-20 person) Plastic Box c/w Bracket - 280 x 280mm	A
34110	HSE - Premium (1-50 person) Plastic Box c/w Bracket - 330 x 340mm	A

Metal Cabinets

Our range of metal cabinets are ideal for first aid supplies storage and are also suitable for clinical environments.

- Supplied with two keys for added security
- Fixed shelves
- Size options available

New glass door metal cabinets available:

- Toughened durable glass door
- Safe and visible supplies
- Quick release push-button
- Secure lock with keys



Content not included



M6 Metal Cabinet - 460 x 300 x 140mm

A

M5 Metal Cabinet - 300 x 300 x 120mm

A



M2 Metal Cabinet - 530 x 530 x 190mm

A



M1 Metal Cabinet - 460 x 300 x 140mm

A

Metal Cabinet First Aid Kits

Robust, lockable and wall mountable metal First Aid Cabinets.

- Available with both HSE and BS8599-1 content
- Supplied with two keys for added security
- Complete with two internal shelves



36406	Metal Cabinet - BS8599-1 Small Content 460 x 300 x 140mm	A
36408	Metal Cabinet - BS8599-1 Medium Content 460 x 300 x 140mm	A
36410	Metal Cabinet - BS8599-1 Large Content 460 x 300 x 140mm	A

70106M	Metal Cabinet - HSE 10 Person 460 x 300 x 140mm	A
70108M	Metal Cabinet - HSE 20 Person 460 x 300 x 140mm	A
70110M	Metal Cabinet - HSE 50 Person 460 x 300 x 140mm	A

Complete First Aid Kit Refills

The HSE Approved Code of Practice states the contents of First Aid containers should be fully restocked after use.

Sufficient supplies should be held in a back-up stock on site. Care should be taken to check the expiry date of each item regularly and replace as required with new stock. Items which are passed the expiry date should be disposed of safely.

- Available with both HSE and BS8599-1 compliant content

32816	HSE 10 Person Refill	A
32818	HSE 20 Person Refill	A
32820	HSE 50 Person Refill	A
32814	HSE Travel Refill	A
32490	BS8599-1 Compliant Small Refill	A
32491	BS8599-1 Compliant Medium Refill	A
32492	BS8599-1 Compliant Large Refill	A

Multipack Replacement Items

36110	Triangular Bandage - 6 pk	A
36102	No.16 Sterile Eye Pad & Dressing - 10 pk	A
36106	Large Sterile Dressing 18 x 18cm - 6 pk	A
36104	Medium Sterile Dressing 12 x 12cm - 6 pk	A
87500	Washproof Assorted Plaster - 120 pk	A
36108	Quick Fix Finger Dressing - 10 pk	A
36204	Microporous Tape 2.5cm x 10m - 12 pk	A
77895	Assorted Safety Pins - 6 pk	A
82124	Moist Cleansing Wipes - 100 pk	A
36130	Disposable Nitrile Gloves - 10 pk	A
36150	Foil Blanket - 10 pk	A
36502	Resuscitator - 10 pk	A
91112	Heavy Duty Clothing Shears	A
36304	Conforming Bandage 7.5cm x 4.5m - 12 pk	A

Vehicle First Aid Kits

Available in either Plastic Box or Nylon Zip pouches for easy storage. These kits contain a comprehensive selection of First Aid products to cover everyday eventualities.

Contents	70140/1	70145/6
First Aid Guidance Leaflet	1	1
Adhesive Plasters Assorted Sizes	20	20
Moist Cleansing Wipes	6	6
Antiseptic Cream	1	1
Microporous Tape	1	1
Conforming Bandage 5cm x 4m	1	1
Conforming Bandage 7.5cm x 4m	1	-
Large Sterile Dressing 18 x 18cm	1	1
Medium Sterile Dressing 12 x 12cm	-	2
Low Adherent Dressing 5 x 5cm	2	2
Low Adherent Dressing 10 x 10cm	-	2
Adhesive Dressing 7cm x 6cm	2	2
Sterile Saline Solution 20ml Pod	-	2
Emergency Foil Blanket	-	1
Safety Pins	6	6
First Aid Scissor	1	1
Triangular Bandage	2	1
Nitrile Gloves	2	1

Alcohol Breath Tests

NF Certified & Approved single use disposable breathalyser kit.

Complies with all necessary requirements to the legislation effective from July 2012 requiring every vehicle in France to carry a minimum of one unused breathalyser at all times.



97084	NF Breath Test Twin Pack	B
-------	-----------------------------	---

Travel and Domestic Kits



70140	Travel and Domestic Kit - Plastic Box 227 x 165 x 47mm	A
70141	Travel and Domestic Kit - Nylon Pouch 180 x 110 x 45mm	A
70145	UK Motoring Kit - Plastic Box 227 x 165 x 47mm	A
70146	UK Motor Kit - Nylon Zip Bag 180 x 110 x 45mm	A

Sterile Medical Equipment Pack

Guarantee UK standards of sterility and protection from risk of infection. Contents:

- 2 Blood Lancets, 23g
- 2 Skin Closure Strips 6mm x 75mm (3)
- 1 Cannula with Port, Green 18g x 45mm
- 2 Non-adherent Dressings 5cm x 5cm
- 2 Non-adherent Dressings 10cm x 10cm
- 20 Plasters, Washproof Assorted
- 1 Microporous Tape 2.5cm x 10m
- 3 Needles, Blue 23g x 1"
- 3 Needles, Green 21g x 1 1/2"
- 3 Needles, Orange 25g x 5/8"
- 4 Pre-injection Swabs
- 1 Certificate
- 1 Giving Set, 185cm
- 1 Gloves (Pair)
- 1 Suture W319
- 2 Syringe 1ml
- 2 Syringe 2.5ml
- 2 Syringe 5ml



70172	Overseas Travellers Pack in Sealed Bag 250 x 140 x 35mm	C
-------	--	---

This range is suitable for carrying in large commercial and light goods service vehicles. Kits include sterile saline solution and some have brackets to allow permanent fixing within the vehicle.

Public Carrying Vehicle First Aid Kits

Contents	70129	70126B	70128
First Aid Guidance Leaflet	1	1	1
First Aid Carried Sticker	-	1	1
No.16 Sterile Eye Pad	1	1	2
No.2 Sterile Ambulance Dressings	3	3	3
Triangular Bandage	2	2	2
Adhesive Plasters Assorted Sizes	30	30	30
Moist Cleansing Wipes	10	10	10
Safety Pins	12	12	12
Disposable Nitrile Glove	1	1	1
First Aid Scissors	1	-	-
Heavy Duty Clothing Shears	-	1	1
Conforming Bandage - 5cm x 4m	-	1	1
Conforming Bandage - 7.5cm x 4m	1	-	1
Sterile Eye Wash Solution - 500ml	-	-	1



70129	PCV Kit - Standard Content Plastic Box 227 x 165 x 47mm	A
70129B	PCV Kit - Standard Content c/w Bracket Plastic Box 227 x 165 x 47mm	A
70128	PCV Kit - Standard Extra c/w Bracket & Eye Wash Plastic Box - 240 x 170 x 80mm	A
70126B	PCV Kit - Standard Extra Content c/w Bracket Plastic Box - 240 x 170 x 80mm	A

Commercial Vehicle First Aid Kits

Supplied with brackets to allow permanent fixing within a vehicle.

Contents	70124B	70130B
First Aid Guidance Leaflet	1	1
No.16 Sterile Eye Pad	-	1
Medium Sterile Dressing 12 x 12cm	-	2
Large Sterile Dressing 18 x 18cm	1	1
Low Adherent Dressing 10 x 10cm	-	2
Triangular Bandage	2	1
Adhesive Plasters Assorted Sizes	-	20
Moist Cleansing Wipes	4	6
Safety Pins	6	12
Disposable Nitrile Glove	1	1
Heavy Duty Clothing Shears	-	1
Sterile Eye Wash Solution Pods - 20ml	2	-
Sterile Eye Wash Solution - 500ml	-	1
Microporous Tape	-	1
Emergency Foil Blanket	1	-

**BEST
SELLER**



70130B	Large Goods Vehicle First Aid Kit Plastic Box c/w Bracket 240 x 170 x 80mm	A
70124B	Light Commercial Vehicle First Aid Kit Plastic Box c/w Bracket 165 x 115 x 55mm	A

INCIDENT & EVACUATION FIRST AID

A range of Incident and Evacuation First Aid Kits designed to manage disaster, mass casualty and extreme trauma situations.

Our kits have a comprehensive content, some with highly engineered products such as Celox and CAT Military tourniquets, for complete assurance adequate cover is available in the event of an emergency.

Trauma First Aid Kit

Grab bag first aid kit for the treatment of major bleeding and burns. The kit includes celox haemostatic granules, which clot blood fast. The kit also includes a selection of high absorption trauma dressings and a quick release tourniquet.



CONTENTS

No.16 Eye Pad	1
Sterile 12 x 12cm Dressing	1
Triangular Bandage	1
Nitrile Disposable Gloves (Pair)	2
Moist Cleansing Wipes	10
Microporous Tape 2.5cm x 10m	1
Sterile Eye Wash Solution 20ml Pods	5
Sterile Gauze Swabs 10 x 10cm (Pack of 5)	1
Trauma Dressing 15 x 18cm	1
Quick Fix Finger Dressing	2
Military Field Dressing 20 x 19cm	1
Burn Dressing 10 x 10cm	1
Quick Release Tourniquet	1
Emergency Foil Blanket	1
Heavy Duty Clothing Shears	1
Celox Haemostatic Granules 35g	1

32944 Trauma First Aid Kit
Grab Bag - 255 x 155 x 95mm

B

Emergency Rucksack First Aid Kit

Our emergency rucksack first aid kit is ideal for administering first aid on the move.



CONTENTS

Guidance leaflet	1
No.16 Sterile Eye Pads	3
Sterile 12 x 12cm Dressings	6
Sterile 18 x 18cm Dressings	2
Washproof Plasters (Pack of 20)	3
Triangular Bandage	3
Nitrile Disposable First Aid Gloves	9
Yellow Disposal Bags	6
Moist Cleansing Wipes	30
Safety Pins (Pack of 6)	2
Microporous Tape 2.5cm x 10m	1
Sterile Eye Wash Solution 20ml Pods	5
Sterile Eye Wash Solution 250ml	2
Conforming Bandage 7.5cm x 4m	2
Sterile QuickFix Finger Dressing	3
Burn Dressing 10 x 10cm	2
Burn Blotts 3.5g	3
Instant Ice Pack	1
Emergency Foil Blanket	2
First Aid Scissor	1
Heavy Duty Clothing Shears	1
Resuscitade	1

32954 Emergency Rucksack First Aid Kit
Rucksack - 380 x 315 x 145mm

B

Remote Medical First Aid Kit

The remote medical kit includes a range of resuscitation devices, burn treatments and splints. The kit includes a cervical collar in cases of suspected head, neck or spinal injuries.



NEW

CONTENTS

Sterile No.16 Eye Pad	2
Triangular Bandage	1
Nitrile Disposable Gloves (Pair)	1
Disposal Bag (Pack 6)	1
Casualty Monitoring Card	2
First Aid Manual	1
Moist Cleansing Wipes	10
Fabric Strapping 2.5cm x 2m	1
Sterile Eye Wash 500ml	1
White Fingerbobs (Pack 6)	1
Blue Visually Detectable Tape 2.5cm x 5m	1
Conforming Bandage 5cm x 4m	1
Conforming Bandage 10cm x 4m	1
Crepe Bandage 5cm x 4.5m	1
Sterile QuickFix Finger Dressing	1
Burn Dressing 10 10cm	1
Burn Dressing Face Mask	1
Cling Film	1
Waterproof Assorted Plasters	100
First Aid Spray 150ml	1
Instant Ice Pack	1
Disposable Towel	2
Emergency Foil Blanket	1
Forceps	1
Heavy Duty Clothing Shears	2
Tweezers	1
Disposable Razors	1
Oxygen Mask Set	1
Single Use Resuscitator	1
Pocket CPR Face Mask	1
Stethoscope	1
Digital Wrist Blood Pressure Monitor	1
Finger Splint	1
Arm Splint	1
Leg Splint	1
Cervical Collar	1
Universal Bite Guard	1

33933

Remote Medical First Aid Kit

First Aid Holdall - 360 x 280 x 180mm

B

Acid Attack Response First Aid Kit

The acid attack response first aid kit contains a range of products for washing chemicals from the skin. Content includes items for treating burns and irrigating eyes following an acid, bleach or chemical attack. Also contains safety equipment protecting the rescuer such as goggles, gauntlets and a face mask.



NEW



CONTENTS

Goggle	1
Class FFP3V Valved Fold - Flat -Singles	1
Nitri-Tech Glove - 13	1
Nitrile Disposable Glove - Pair	6
Bottle of Water - 500ml	6
Litmus Test Strips - Pack of 100	1
Eye Wash Solution 20ml - Singles	10
Eye Wash Solution 500ml	1
Gauze Swabs - 7.5 x 7.5cm (100)	1
Burn Dressing 20 x 20cm	2
Emergency Foil Blanket - 2.1m x 1.3m	3
Heavy Duty First Aid Snips	1

32992

Acid Attack Response First Aid Kit

Active Holdall Bag Orange - 380 x 260 x 170mm

B

INCIDENT & EVACUATION FIRST AID

Our incident first aid kits have a comprehensive content in case of a serious incident. The kit provides products for the treatment of serious bleeds, eye injuries and burns.

CONTENTS	33904	32966 / 32970	32997
Guidance Leaflet	-	-	1
Triangular Bandage	6	4	8
Sterile No.16 Eye Pad	-	-	6
Sterile Medium First Aid Dressing 12 x 12cm	-	-	12
Sterile Large First Aid Dressing 18 x 18cm	-	2	4
Conforming Bandage 10cm x 4m	-	-	6
Crepe Bandage 10cm x 4.5m	-	3	-
Nitrile Gloves (Pair)	20	3	10
Microporous Tape 2.5cm x 10m	10	1	1
Sterile Eye Wash Solution 20ml Pods	10	10	25
Sterile Eye Wash Solution 500ml	2	-	-
Moist Cleansing Wipes	100	10	100
Sterile Ambulance Dressing No.3 28 x 18cm	12	-	4
Sterile Ambulance Dressing No.4 33 x 20cm	-	-	4
Trauma Dressing 10 x 18cm	5	-	-
Trauma Dressing 15 x 18cm	5	-	-
Military Field Dressings	-	3	-
Sterile Low Adherent Dressing 5 x 5cm	-	-	5
Sterile Low Adherent Dressing 7.5 x 7.5cm	-	-	5
Sterile Low Adherent Dressing 10 x 10cm	-	-	5
Tearable Zinc Oxide Tape 2.5cm x 10m	4	-	-
Emergency Foil Blanket	10	2	6
Scissors Blunt/Sharp	2	-	-
Heavy Duty Clothing Shears	-	1	1
Resuscitator	10	1	10
Clinell Hand Gel 250ml	1	-	-
Glucogel Hypostop Pk 3	-	2	-
Combat Application Tourniquet	-	1	-
Quick Release Tourniquet	5	-	-
Celox Granules 15g	-	-	1
Celox FoxSeal ChestSeal Pk2	-	-	1
Premium Burns Medium Kit	-	-	1
BurnSoothe Dressing 10 x 10cm	-	2	-
Safety Pins	-	-	12
Wound Closure Strips 3 x 75mm (Strip of 5)	-	-	6
Assorted Fabric Plasters Box 120	-	-	1
Washproof Plasters	-	30	-
All weather notebook	-	1	-
HB Rubber Tipped Pencil	-	1	-
Instant Ice Pack	-	2	-
Flexi Splint Leg - 11x100cm	-	1	-
Gauze Swabs Sterile - 10cm x 10cm (5)	-	2	-

32970	Incident First Aid Kit Active Holdall Bag Orange - 380 x 260 x 170mm	B
32966	Incident First Aid Kit ABS Box & Wall Bracket - 410 x 300 x 150mm	B
33904	Critical Incident First Aid Kit First Aid Holdall - 360 x 280 x 180mm	B
32997	Major Incident First Aid Kit Trauma Bag - 410 x 350 x 285mm	B

32970



32966



33904



32997



Mass Casualty First Aid Kit

Our mass casualty first aid kit provides access to specialist provisions in the event of a major incident. The kit provides provisions for a range of injuries including burns, lacerations, eye injuries and respiratory problems.



CONTENTS

Sterile No.16 Eye Pads	10
Nitrile Disposable Gloves (Pair)	10
Microporous Tape 2.5cm x 10m	10
Low Adherent Dressings 10 x 10cm (Box 25)	1
Conforming Bandage 7.5cm x 4m	10
Sterile Gauze Swabs 10 x 10cm (5)	5
Ambulance Dressing No.2 (20 x 17.5cm)	5
Ambulance Dressing No.3 (28 x 18cm)	5
Ambulance Dressing No.4 (33 x 20cm)	5
Finger Dressing	10
Military Field Dressing 10 x 19cm	4
Military Field Dressing 20 x 19cm	1
Burn Dressing 10 x 10cm	5
Resusciate	10
Clinell Hand Gel 100ml	1
CAT Military Tourniquet	1
Celox Rapid 5ft Z-Fold	1
Sterile Eye Wash 20ml Pods	10

32999	Mass Casualty First Aid Kit Trauma Bag - 410 x 350 x 285mm	B
-------	---	---

Customisation

We offer customised first aid kits should you have a specific requirement or content in mind. Should you require further information about any of the products please contact our sales office.

CELOX

Celox is a range of haemostats which control moderate to severe traumatic bleeding.

Proven technology that is applied in military medicine and by the emergency services. Choose from our range of Celox haemostatic granules and Celox gauze dressings.



99806	Celox Granules 15g	B
99803	Celox Gauze 5ft Z-Fold	B
99804	Celox Rapid 5ft Z-Fold	B
99811	Celox FoxSeal Chest Seal - Pk 2	B

REDCAP

Redcap™ chemical eye wash solution 500ml is a phosphate buffer solution and has a neutralising effect.

Redcap™ can be used to neutralise both acids and alkali's on contact returning eyes and skin to its natural pH value.



80270	Chemical Eye Wash Solution 500ml c/w Eye Bath	B
80272	Redcap™ Chemical Splash Eye Wash Station	B
80274	Redcap™ Chemical Eye Wash Solution 500ml	B

Forestry First Aid Kits

HSE Workplace First Aid guidelines apply to any workplace, both indoors and out.

Employers are required to provide comprehensive First Aid cover for workers in all situations and should be available to employees at all times.

Active Forestry First Aid Kit

The Active Forestry First Aid Kit has a content suitable for a smaller team of workers. The Active bag is compact and versatile; can be carried via shoulder or belt strap when working remotely. The belt strap folds neatly away behind a Velcro pad when not required.

Durable, rip resistant and easy to clean our Active Bag is extremely practical in all pursuits.

- Internal, external side pockets
- Adjustable shoulder and belt straps
- Secure key holder
- Strong, water-resistant, highly visible material



70165	Active Forestry First Aid Kit Active Bag - Green 300 x 175 x 100mm	A
701650	Active Forestry First Aid Kit Active Bag - Orange 300 x 175 x 100mm	A

Contents	70160	70165/O	70168/O
First Aid Guidance Leaflet	1	1	1
Adhesive Plasters Assorted Sizes	1	2	2
Non Woven Triangular Bandage	1	2	6
Nitrile Disposable Gloves	1	3	6
Sterile Moist Cleansing Wipes	2	6	10
Assorted Safety Pins (pk6)	1	1	1
Ambulance Dressing No.4	1	1	2
Ambulance Dressing No.3	-	1	-
Trauma Dressing 15 x 18cm	-	-	1
Trauma Dressing 10 x 18cm	-	1	-
Resusciate	1	2	1
Sterile Eye Dressing No.16	-	2	2
Sterile Medium Dressing 12x12cm	-	2	6
Sterile Large Dressing 18x18cm	-	1	4
Sterile No.3 X-Large Dressing 28x17cm	-	-	3
Sterile Eye Wash Solution 20ml Pods	-	5	5
Sterile Eye Wash Solution 500ml Bottles	-	-	1
Conforming Bandage 7.5cmx4m	-	-	1
Foil Blanket	-	1	1
Heavy Duty First Aid Snips	-	1	1
Plastic Tweezers	-	1	1

Personal Forestry First Aid Kits

Current Forestry Commission guidelines require a Forestry and Chainsaw worker to carry a personal First Aid Kit which includes a large wound dressing.

This kit comes in a belt pouch which can easily strap to a belt or harness.



70160	Personal Forestry First Aid Kit Nylon Pouch - 180 x 110 x 45mm	A
--------------	--	----------

Forestry & Chainsaw First Aid Kits

MediKit's Forestry and Chainsaw First Aid Kits provide comprehensive content to cover a troop of workers with more dressings, gloves and eye wash available than a standard First Aid kit.

Available in a standard first aid holdall or alternatively in the Active range first aid holdall complete with protective padded lining and high visibility, rip resistant material.



70168	Forestry and Chainsaw First Aid Kit Large Green Holdall 360 x 280 x 180mm	A
701680	Forestry and Chainsaw First Aid Kit Active Large Orange Holdall 380 x 260 x 170mm	A
70168G	Forestry and Chainsaw First Aid Kit Active Large Green Holdall 380 x 260 x 170mm	A

Farming and First Response Kits

Tractor First Aid Kits

Supplied with brackets to allow permanent fixing within a vehicle.

Contents	70124	32125
First Aid Guidance Leaflet	1	1
Non Woven Triangular Bandage	2	2
Nitrile Disposable Gloves	1	1
Sterile Moist Cleansing Wipes	4	4
Assorted Safety Pins (pk6)	1	1
Sterile Large Dressing 18x18cm	1	1
Sterile Eye Wash Solution 20ml Pods	2	2
Foil Blanket	1	1



70124	Tractor First Aid Kit Plastic Box - 165 x 115 x 55mm	A
32125	Tractor First Aid Kit Nylon Pouch - 180 x 110 x 45mm	A

Contents	70918	70153FB
First Aid Guidance Leaflet	1	1
No.16 Sterile Eye Pad	2	6
Medium Sterile Dressing 12 x 12cm	6	12
Large Sterile Dressing 18 x 18cm	2	4
Adhesive Plasters Assorted Sizes	40	60*
Triangular Bandage	6	6
Disposable Nitrile Gloves	6	6
Moist Cleansing Wipes	20	30
Safety Pins	12	12
Microporous Tape 2.5cm x 10m	1	-
Sterile Eye Wash Solution 20ml Pods	5	2
Conforming Bandage 7.5cm x 4m	1	1
Crepe Bandage 7.5cm x 4.5m	2	-
Low Adherent Dressing 10 x 10cm	4	-
Sterile Finger Bandage	2	2
Burn Dressing 10 x 10cm	1	1
Hydrogel Burn Blotts 10ml	-	3
Instant Ice Pack	2	-
Emergency Foil Blanket	3	1
Heavy Duty Clothing Shears	1	1
Resuscitator Face Shield	1	3
Blue Fingerstalls	-	6
Blue Detectable Washproof Tape	-	1

*Blue Detectable Plasters Supplied in Food Processing / Catering Kits

First Response Kit

First Response First Aid Kits are compliant with the BS 8599-1 Standard. They have additional items suitable for treating patients who may be some distance from a static First Aid Station.

- Shoulder strap
- Multiple compartments for easy storage
- The bag has extra space which can be used to include additional products



70918	First Response Kit Soft Nylon Shoulder Bag - 280 x 230 x 120mm	A
--------------	--	----------

Catering First Response Kit

BS 8599-1 compliant plus additional items including blue adhesive plasters and burn blotts.

- Robust ABS plastic boxes with internal compartments
- Rubber lid seal preventing contamination from dust and moisture
- Complete with clip in wall bracket



70153FB	Catering First Response Kit ABS Plastic Box - 320 x 220 x 130mm	A
----------------	---	----------

First Aid for Burns

MediKit's range of general purpose Burns Kits are suitable for small working environments; each kit includes various sizes of hydrogel burn dressings, along with standard first aid items. Choose from a selection of boxes and bags.

Burns Kits

84050	Premium Burns Kit Medium Content - 280 x 280mm	B
84052	Premium Burns Kit Large Content - 280 x 280mm	B
34846	Standard Burns Kit Medium Content - 280 x 280mm	B
34847	Standard Burns Kit Large Content - 280 x 280mm	B



Contents	Medium	Large	Responder	Deluxe Responder
First Aid Guidance Leaflet	1	1	-	1
BS8599-1 Small Refill	-	-	1	-
Medium Sterile Dressing 12 x 12cm	1	2	-	-
Large Sterile Dressing 18 x 18cm	1	2	-	-
Disposable Nitrile Gloves	2	4	-	2
Safety Pins 6	1	2	-	-
Burn Dressing 10 x 10cm	1	1	1	8
Burn Dressing 20 x 20cm	1	2	1	3
Burn Dressing 60 x 40cm	-	-	-	2
Burn Limb Dressing 10 x 40cm	-	-	-	2
Burn Blotts 10ml	2	4	2	3
Burn Gel 50ml	1	1	1	-
Burn Gel 125ml	-	-	-	3
Heavy Duty Clothing Shears	1	1	-	1
Underwrap Foam Bandage 7cm x 27m	-	-	-	2
Emergency Foil Blanket	-	-	-	1
Conforming Bandage 7.5cm x 4m	-	-	-	2
Microporous Tape - 2.5cm x 10m	-	-	-	1

Paramedic Bum Bag Burns Kits

The Paramedic Bum Bags fold out to reveal multiple compartments for easy organisation of supplies. The material is durable and wipe clean, making it practical for use in many working environments.

34846PB	Paramedic Bum Bag Burns Kit Medium Content - 270 x 160 x 45mm	B
34847PB	Paramedic Bum Bag Burns Kit Large Content - 270 x 160 x 45mm	B



First Aid for Burns

Burnshield Products

Burnshield is a simple to use gel impregnated burn treatment designed to lower skin temperature and alleviate pain prior to hospital treatment.

- Available as dressings or blankets. Suitable for treating first, second and third degree burns and scalds.
- Manufactured from open cell foam, all dressings come saturated in a water based gel made sterile by gamma irradiation.
- Simple to apply and remove without causing any discomfort, Burnshield soothes and moisturises the burn area by inhibiting fluid loss.
- All items may be purchased individually to replenish used stock.



84070	Burn Blott 3.5ml - Pack of 10	B
84030	Burnshield - 10cm x 10cm	B
84032	Burnshield - 20cm x 20cm	B
84034	Burnshield - 60cm x 40cm	B
84040	Burnshield Blanket	C
84074	Burnshield Hydrogel - 50ml	B
84075	Burnshield Hydrogel -125ml	B

Burnshield Hydrogel

- Absorbs and dissipates heat - minimises burn damage
- Cooling - pain relief and shock reduction
- Clear - easy to assess wound
- Non-stick - will not adhere to wound
- Natural anti-bacterial
- Non Toxic
- Non Irr

Burnshield Dressings

- Foam based dressing for optimal heat dissipation.
- Conforms to wound area
- Fibreless elastic opencell foam sheeting
- Wound covering

Burns Responder Kits

First Response First Aid Kits comply with the BS8599-1 Standard plus many additional items which could be required when treating a patient who may be some distance from a static First Aid Station.

MediKit's Deluxe Burns Responder kit has the greatest variety of hydrogel dressings in different sizes, all supplied in our first aid holdall for greater mobility.

- Detachable shoulder straps
- Internal compartments for easy storage and access.
- The bags have extra space which can be used to include additional products.

84056	Burns Responder Kit Soft Nylon Shoulder Bag 280 x 230 x 120mm	B
84345	Deluxe Burns Responder Kit Large Green Holdall 360 x 280 x 180mm	B



Sports First Aid Kits

MediKit's Sports First Aid Kits are supplied in various content sizes with different styles of bags to suit multiple applications. We are able to customise kits to meet individual requirements, please contact our Sales Office for more information on all the options available.

70190	Sports Player First Aid Kit Nylon Zip Pouch - 210 x 140 x 60mm	B
70189	Sports Team First Aid Kit Active Bag Green - 175 x 300 x 100mm	B
701890	Sports Club First Aid Kit Active Bag Orange - 175 x 300 x 100mm	B
70188	Sports Club First Aid Kit Active Holdall Bag Green - 380 x 260 x 170mm	B
701880	Sports Club First Aid Kit Active Holdall Bag Orange - 380 x 260 x 170mm	B

Contents	70190	70189/0	70188/0
Guidance Card	1	1	1
Adhesive Plasters Assorted Sizes	20	20	40
Blister Plaster Pack	1	1	1
Safety Pins	6	6	6
Moist Cleansing Wipes	4	10	20
No.16 Sterile Eye Pad	-	1	1
Medium Sterile Dressing 12 x 12cm	-	1	1
Large Sterile Dressing 18 x 18cm	-	1	2
Sterile Swabs 5 x 5cm	-	-	5
Trauma Dressing 10 x 18cm	-	1	1
Trauma Dressing 15 x 18cm	-	-	1
Nitrile Gloves	2	2	5
Microporous Tape	1	1	1
Zinc Oxide Tape 5cm x 10m	-	-	1
Tearable Zinc Oxide Tape 5cm x 10m	-	-	1
Antiseptic Cream 30g	1	1	1
Triangular Bandage	-	2	2
Conforming Bandage 5cm x 4m	1	-	1
Conforming Bandage 7.5cm x 4m	-	1	1
Crepe Bandage 5cm x 4.5m	-	-	1
Crepe Bandage 7.5cm x 4.5m	-	1	1
Cohesive Bandage 5cm x 4.5m	-	1	2
Cohesive Bandage 7.5cm x 4.5m	-	-	1
Low Adherent Dressing 5 x 5cm	2	1	3
Low Adherent Dressing 7.5 x 7.5cm	-	1	2
Adhesive Dressing 7 x 6cm	-	1	2
Adhesive Dressing 8.3 x 6cm	2	1	2
Adhesive Dressing Strip 7.5cm x 1m	1	-	1
EAB 5cm x 4.5m	-	-	1
Tearable EAB 5cm x 4.5m	1	1	1
Underwrap Roll	-	1	1
Instant Ice Pack	-	2	4
Scissor	1	-	-
Heavy Duty Clothing Shears	-	1	1
First Aid Spray - 70ml	1	1	-
First Aid Spray - 400ml	-	-	1
Resusciate	1	1	2
Wound Closure - Strip of 5	1	1	2
PPP Pack	-	-	1
Plastic Tweezers	-	1	1
Foil Blanket	-	1	2
Eye Wash Pods - 20ml	-	2	2
Petroleum Jelly - 190g	-	-	1
Muscle Rub	-	-	1

Player First Aid Kit

Basic sports kit content contained in an easily accessed first aid pouch with belt loop.



Team First Aid Kit

Medium sized sports kit containing a wide range of sports consumables, supplied in a versatile first aid bag.

- Compact bag can be used as a shoulder bag, bum bag or simply carried by hand
- Durable, rip resistant and highly visible



Club First Aid Kit

Large, comprehensive content containing our widest selection of sports consumables including muscle rub, tearable tapes and trauma dressings for more serious injuries.

The Club First Aid Kit is supplied in a first aid holdall with multiple compartments and pockets for easy organisation of supplies.

- Durable, rip resistant and highly visible
- Protective padded lining
- Detachable internal compartments



General Sport First Aid Kits

General Sports First Aid Kits

A basic range of First Aid Kits suitable for general sports injuries.



Contents	70192	70194/7	70198
Guidance Card	1	1	1
Adhesive Plasters Assorted Sizes	20	40	120
Safety Pins	6	6	6
Moist Cleansing Wipes	10	20	20
No. 16 Sterile Eye Pad	-	1	1
Medium Sterile Dressing 12 x 12cm	1	1	1
Large Sterile Dressing 18 x 18cm	1	2	2
Nitrile Gloves	2	2	10
Microporous Tape	1	1	1
Antiseptic Cream 30g	1	1	1
Triangular Bandage	2	2	2
Conforming Bandage 5cm x 4m	1	2	-
Conforming Bandage 7.5cm x 4m	1	1	2
Conforming Bandage 10cm x 4m	-	-	2
Crepe Bandage 5cm x 4.5m	-	1	2
Crepe Bandage 7.5cm x 4.5m	-	1	2
Low Adherent Dressing 5 x 5cm	2	3	-
Low Adherent Dressing 7.5 x 7.5cm	1	2	25
Adhesive Dressing 7 x 6cm	-	-	2
Adhesive Dressing 8.3 x 6cm	2	1	1
Adhesive Fabric Dressing Strip 7.5cm x 1m	1	1	1
EAB 5cm x 4.5m	1	1	1
EAB 7.5cm x 4.5m	-	-	2
Instant Ice Pack	1	3	10
Scissor	1	1	1
Heavy Duty Clothing Shears	-	-	1
First Aid Spray - 70ml	1	-	-
First Aid Spray - 150ml	-	1	1
Resusciate	1	2	2
Wound Closure - Strip of 5	2	2	2
PPP Pack	-	1	1
Foil Blanket	-	-	6
Eye Wash Pods - 20ml	-	-	25
Petroleum Jelly - 250ml	-	-	1

Delph Sun Protection

Delph Sun Care offers a range of high quality products from SPF 15 to SPF 50 with a 4* superior UVA and UVB protection. Delph products are fully dermatologically tested to medical standards (on humans) and are suitable for all skin types.

90164	Delph Sun Lotion SPF30 c/w Pump - 2 ltr	D
90170	Delph Sun Lotion SPF15 - 150ml	C
90171	Delph Sun Lotion SPF30 - 150ml	C



70192	Sports Kit - Medium Nylon Zip Pouch 210 x 140 x 60mm	B
70194	Sports Kit - Large Grab Bag 280 x 230 x 120mm	B
70197	Sports Kit - Large Sports Holdall 360 x 280 x 180mm	B
70198	Sports Kit - Comprehensive Sports Holdall 360 x 280 x 180mm	B

Petroleum Jelly

90502	Petroleum Jelly - 190g	B
90504	Petroleum Jelly - 400g	B
90506	Petroleum Jelly - 2kg	B
90508	Petroleum Jelly - 5kg	B



90172	Delph Sun Lotion SPF50 - 150ml	C
90174	Delph After Sun Lotion - 150ml	C
90175	Delph After Sun Gel - 150ml	C

Sport First Aid Consumables

Zinc Oxide Tape


- Non-stretch cotton fabric
- Powerful, long lasting adhesive for immobilising joints
- Designed to tear by hand



89008	Zinc Oxide Tape - 2.5cm x 10m	C
89009	Zinc Oxide Tape - 3.8cm x 10m	C
89010	Zinc Oxide Tape - 5cm x 10m	C
89004	Tearable Zinc Oxide Tape - 2.5cm x 10m	C
81652	Tearable Zinc Oxide Tape - 3.8cm x 10m	C
81654	Tearable Zinc Oxide Tape - 5cm x 10m	C

EAB - Non Tearable Tape

- High strength adhesive
- Strongly woven for firm support



80120	Elastic Adhesive Bandage - 2.5cm x 4.5m	B
80122	Elastic Adhesive Bandage - 5cm x 4.5m	B
80124	Elastic Adhesive Bandage - 7.5cm x 4.5m	B
80126	Elastic Adhesive Bandage - 10cm x 4.5m	B

EAB - Tearable Tape



89024	Elastic Adhesive Bandage Tearable - 5cm x 4.5m	C
89026	Elastic Adhesive Bandage Tearable - 7.5cm x 4.5m	C

Cohesive Bandages

- Latex free, durable and strong
- Sticks to itself, not skin or hair
- Tearable by hand



80142	Cohesive Bandages - 5cm x 4.5m	C
80144	Cohesive Bandages - 7.5cm x 4.5m	C
80145	Cohesive Bandages - 10cm x 4.5m	C


Foam Bandage

- Extra comfort and protection
- Applied prior to taping to prevent chafing
- Thin, lightweight foam bandage



89052	Underwrap foam bandage - 7cm x 27m	C
-------	------------------------------------	---

Dressing Strip



89062	Fabric Adhesive Dressing Strip - 7.5cm x 5m	C
89054	Fabric Adhesive Dressing Strip - 6cm x 5m	C
81610	Fabric Adhesive Dressing Strip - 7.5cm x 1m	C

Kinesiology Tape

Specialised therapeutic tape designed to help you overcome injury and pains; it works with the body to support recovery without limiting movement.

- High stretch and flexible tape
- Breathable and hypo-allergenic wave adhesive



81626	Kinesiology Tape Light Tan - 5cm x 5m	C
81627	Kinesiology Tape Light Blue - 5cm x 5m	C
81628	Kinesiology Tape Black - 5cm x 5m	C

Treatments



94058	Supercool Freeze Gel - 100ml tube	C
90657	Vapour Rub - 50g	C
90417	Massage Oil - 500ml	C
90659	Muscle Rub - 200g	C
90663	Sports Balm - 500ml	C

Physio Felt

89002	Adhesive Physio Felt 5mm - 22.5 x 45cm	C
89003	Adhesive Physio Felt 7mm - 22.5 x 45cm	C

Hot And Cold Treatment



90180	Instant Cold Pack - Small	C
90182	Instant Cold Pack - Standard	C
90188	Reusable Hot / Cold Pack	C
90116	First Aid Cool Spray - 150ml	C
90118	First Aid Cool Spray - 400ml	C
90119	First Aid Heat Spray - 150ml	C

Miscellaneous Tape

General purpose sports tapes.



81641	Sport PVC Sock Tape - 1.9cm x 20m	C
81642	Dressing Retention Tape - 5cm x 10m	C
80112	Microporous Tape - 1.25cm x 10m	A
80114	Microporous Tape - 2.5cm x 10m	A

Plasters

Washproof Plasters



87141	4cm x 4cm - Box of 50	B
87143	7.5cm x 2.5cm - Box of 100	B
87145	7.5cm x 5cm - Box of 50	B
87157	Anchor - Box of 50	B
87914	Extension - Box of 50	B
87910	Fingertip - Box of 50	B
87150	Assorted - Box of 100 (5 Sizes)	B

Stretch Fabric Plasters



86221	4cm x 4cm - Box of 50 B	C
86225	7.5cm x 5cm - Box of 50	B
86223	7.5cm x 2.5cm - Box of 100	B
86228	Assorted - Box of 100 (5 Sizes)	B
86230	Extension - Box of 50	B

Premium Fabric Plasters

Premium Fabric plasters are recommended where maximum protection and a long wearing, durable plaster is required.

86021	4cm x 4cm - Box of 100	B
86023	7.5cm x 2.5cm - Box of 100	B
86025	7.5cm x 5cm - Box of 50	B
86038	Assorted - Box of 100 (5 Sizes)	B
86231	Fingertip - Box of 50	B
86232	Anchor - Box of 50	B

Plasters - Assorted Box 120

This quality range is available as a box of 120 assorted plasters in 7 different sizes. They are supplied in durable polypropylene translucent boxes with secure fastenings to protect the content and also to reduce the risk of contamination by dust and moisture.

- Hypoallergenic
- Individually Wrapped and Sterile
- Highly Absorbent Low Adherent Pad



85500	MediKit Blue Detectable Plasters Assorted - Box of 120	A
86500	MediKit Fabric Plasters Assorted - Box of 120	A
87500	MediKit Washproof Plasters Assorted - Box of 120	A

Blue Metal Detectable Plasters



85601	4cm x 4cm - Box of 50	B
85603	7.5cm x 2.5cm - Box of 100	B
85605	7.5cm x 5cm - Box of 50	B
85913	Anchor - Box of 50	B
85910	Fingertip - Box of 50	B
85914	Extension - Box of 50	B
85610	Assorted - Box of 100 (5 Sizes)	B
85636	Premium Assorted - Box of 100 (5 Sizes)	B

Coverplast Adhesive Dressings

- Bacterial and Viral barrier
- Low allergy acrylic adhesive
- Breathable polyurethane film
- Impermeable to oil and fat
- Good skin compatibility
- Highly conformable
- Wide all-round adhesive margin




Coverplast Barrier Washproof

81016	Coverplast Barrier 3.8cm x 3.8cm - Box of 100	C
81022	Coverplast Barrier 6.3cm x 2.2cm - Box of 100	C
81018	Coverplast Barrier 6.3cm x 3.8cm - Box of 100	C
81014	Coverplast Barrier 7.2cm x 2.2cm - Box of 100	C
81020	Coverplast Barrier 7.2cm x 5cm - Box of 100	C
81024	Coverplast Barrier Assorted Sizes - Box of 100	C

Coverplast Classic Washproof

81374	Coverplast Classic 3.8cm x 3.8cm - Box of 100	C
81364	Coverplast Classic 5cm x 4.4cm (Fingertip) - Box of 50	C
81372	Coverplast Classic 7.2cm x 2.2cm - Box of 100	C
81362	Coverplast Classic 7.2cm x 3.8cm - Box of 100	C
81370	Coverplast Classic 7.2cm x 5cm - Box of 100	C
81378	Coverplast Classic Assorted Sizes - Box of 126	C



Anchor




Fingertip



7.5cm x 2.5cm



4cm x 4cm



7.5cm x 5cm
(not to scale)



6.3cm x 3.8cm

Coverplast Barrier Washproof

- Visual Inspection: Signal Blue Contrast Colour
- Metal Detection: Iron Filaments
- Metal Detection: Aluminium Filaments
- X-Ray Detection: Tungsten Filaments

81244	Coverplast Detectable 3.8cm x 3.8cm - Box of 100	C
81242	Coverplast Detectable 7.2cm x 2.2cm - Box of 100	C
81246	Coverplast Detectable 7.2cm x 5cm - Box of 100	C
81248	Coverplast Detectable 5cm x 4.4cm Fingertip - Box of 50	C
81240	Coverplast Detectable Assorted Sizes - Box of 95	C
77767	Coverplast Detectable Blue Plasters (Wallet Pack Of 21)	C

Bandages

Our comprehensive range of Bandages provide compression and support to minor injuries.

These products can be added to your First Aid Kit making the kit more comprehensive.

Limb Size	Child Limbs	Adult Arms	Adult Ankles	Adult Legs	Adult Thighs	Adult Body
Small	A	B	C	D	E	J
Medium	B	C	D	E	F	K
Large	C	D	E	F	G	L



Tubular Bandage Elasticated

- Supportive, comfortable, radial bandage
- Ideal for soft tissue injuries, dislocations and sprains
- Range of sizes for all parts of the body
- Ten metre sizes supplied in dispenser box



81700	Tubular Bandage Size B 6.25cm x 1m Roll	B
81702	Tubular Bandage Size C 6.75cm x 1m Roll	B
81704	Tubular Bandage Size D 7.5cm x 1m Roll	B
81706	Tubular Bandage Size E 8.75cm x 1m Roll	B
81708	Tubular Bandage Size F 10cm x 1m Roll	B
81709	Tubular Bandage Size A 4.5cm x 10m Roll	B
81710	Tubular Bandage Size B 6.25cm x 10m Roll	B
81712	Tubular Bandage Size C 6.75cm x 10m Roll	B
81714	Tubular Bandage Size D 7.5cm x 10m Roll	B
81716	Tubular Bandage Size E 8.75cm x 10m Roll	B
81718	Tubular Bandage Size F 10cm x 10m Roll	B
81719	Tubular Bandage Size G 12cm x 10m Roll	B

Tubular Gauze Bandage

- Supportive, comfortable, radial bandage
- Ideal for soft tissue injuries, dislocations and sprains
- Range of sizes for all parts of the body

81720	Tubular Gauze Bandage Size 01 - 20m Roll	B
81722	Tubular Gauze Bandage Size 12 - 20m Roll	B
81724	Tubular Gauze Bandage Size 34 - 20m Roll	B
81726	Tubular Gauze Bandage Size 56 - 20m Roll	B
81728	Tubular Gauze Bandage Size 78 - 20m Roll	B

Applicators For Bandages

81732	Applicators - Small - Size 01	B
81734	Applicators - Medium - Size 12	B
81736	Applicators - Large - Size 34	B
81738	Applicators - XL - Size 56/78	B

Gauze Swabs

- Individually Packed
- Available Sterile or Non Sterile

83102	Gauze Swabs - 5 x 5cm (100)	B
83104	Gauze Swabs - 7.5 x 7.5cm (100)	B
83106	Gauze Swabs - 10 x 10cm (100)	B
83112	Sterile Gauze Swabs - 5 x 5cm (5)	B
83114	Sterile Gauze Swabs - 7.5 x 7.5cm (5)	B
83116	Sterile Gauze Swabs - 10 x 10cm (5)	B

Conforming Bandage

- Soft, lightweight, conformable
- Superb stretch, comfortable to wear
- Excellent for securing dressings

81414	Conforming Bandage - 5cm x 4m	B
81416	Conforming Bandage - 7.5cm x 4m	B
81418	Conforming Bandage - 10cm x 4m	B
81420	Conforming Bandage - 15cm x 4m	B

Cohesive Bandages

- Non-slip offering excellent support and compression
- Will not stick to skin or hair
- Underwrap not required and tearable by hand
- Individually wrapped

80141	Cohesive Bandage - 2.5cm x 4.5m	B
80142	Cohesive Bandage - 5cm x 4.5m	B
80144	Cohesive Bandage - 7.5cm x 4.5m	B
80145	Cohesive Bandage - 10cm x 4.5m	B

Elastic Adhesive Bandage

- High strength adhesive and heavy weight fabric
- Strongly woven for firm support
- Excellent for weak joints, sprains etc



80120	Elastic Adhesive Bandage - 2.5cm x 4.5m	B
80122	Elastic Adhesive Bandage - 5cm x 4.5m	B
80124	Elastic Adhesive Bandage - 7.5cm x 4.5m	B
80126	Elastic Adhesive Bandage - 10cm x 4.5m	B

Crepe Bandage

- Strong weave and high elasticity
- Superb stretch, comfortable to wear
- Excellent for support, sprains etc

81422	Crepe Bandage - 5cm x 4m	B
81426	Crepe Bandage - 7.5cm x 4m	B
81430	Crepe Bandage - 10cm x 4m	B
81432	Crepe Bandage - 15cm x 4m	B

Blue Visually Detectable First Aid Items



81256	Blue Visually Detectable Adhesive Tape 2.5cm x 5m	B
81078	Blue Finger Bobs - Box of 100	B
80180	Blue Finger Stalls XL - Pack of 10	B
80182	Blue Finger Stalls Large - Pack of 10	B
80183	Blue Finger Stalls Medium - Pack of 10	B
80178	White Finger Bobs - Box 100	B

Cotton Wool

- Various sizes available

83003	Cotton Wool - 100gm	B
83006	Cotton Wool - 500gm	B
83024	Cotton Buds (200)	B
83025	Johnson's Cotton Buds (200)	B

Dressings and Splints

First Aid Dressings are provided to assist the blood clotting process and to stop any foreign bodies from entering the wound which could cause an infection.

These products can be added to your First Aid Kit making the kit more comprehensive.

Low Adherent Dressing Sterile

- Perforated surface prevents the pad sticking to the wound
- Can be used for wounds, minor lacerations or grazes
- Easily removed and highly absorbent



81050	Low Adherent Dressings 5 x 5cm - Box of 25	A
81054	Low Adherent Dressings 7.5 x 7.5cm - Box of 25	A
81055	Low Adherent Dressings 10 x 10cm - Box of 25	A

Triangular Bandages

- Ideal for slings and splinting
- Calico can be washed and reused as appropriate

32896	Non Woven Triangular Bandage 127 x 90cm	A
77737	Sterile Triangular Bandage 127 x 90cm	B

Ambulance Dressings

- Sterile and individually wrapped
- Thicker than standard wound pads for great absorption
- Conforming bandage

83226	Ambulance Dressing No.1 Sterile 12.5 x 10cm	B
83228	Ambulance Dressing No.2 Sterile 20 x 17.5cm	B
83230	Ambulance Dressing No.3 Sterile 28 x 18cm	B
83232	Ambulance Dressing No.4 Sterile 33 x 20cm	B

Large Adhesive Dressing Sterile

- Soft and pliable adhesive dressing
- Can be used for wounds, minor lacerations or grazes
- Highly absorbent, low adherent island wound pad



81090	Large Adhesive Wound Dressings 6 x 7cm - Box of 25	A
81092	Large Adhesive Wound Dressings 8.6 x 6cm - Box of 25	A
81095	Large Adhesive Wound Dressings 10 x 9cm - Box of 25	A
81100	Large Adhesive Wound Dressings 15 x 9cm - Box of 25	B

Sterile Unmedicated Dressing

Sterile dressing with conforming stretch bandage are used as the standard first aid dressing

- Available in a range of sizes



32504	Medium Sterile Dressing - 12 x 12cm	A
32506	Large Sterile Dressing - 18 x 18cm	A
32502	No.16 Sterile Eye Pad	A
77522	No.16 Boxed Sterile Eye Pad	A
83244	Quick Fix Sterile Finger Dressing	B
77880	No.7 Small Sterile Dressing - 5 x 5cm	B
77510	No.8 Medium Sterile Dressing - 10 x 8cm	B
77512	No.9 Large Sterile Dressing - 13 x 9cm	B
77514	No.3 XLarge Sterile Dressing - 28 x 17cm	B

Military Field Dressing

- Vacuum compressed packaging
- Sterile single use dressing
- Advanced fluid absorbency and protection
- Camo packaging



83304	Military Field Dressing - 10 x 19cm	B
83306	Military Field Dressing - 20 x 19cm	B

Sterile Skin Closures

- Skin suture to help secure more severe cuts
- Sealed in sterile transparent envelopes
- Hypo-allergenic adhesive and breathable material



80102	Skin Closure - 3 x 75mm - Pk 50	C
80103	Skin Closure - 6 x 75mm - Pk 30	C
80106	Steristrip Skin Closure - 3 x 75mm - Pk 60	C

Trauma Dressing

- Highly absorbent, non-adherent wound pad
- High pressure bandage
- Super-fast application, quick-fix fastening
- Sterile single use dressing



83236	Trauma Dressing - 10 x 18cm	B
83238	Trauma Dressing - 15 x 18cm	B

Splint

- Immobilise broken bones and sprain injuries
- Lightweight, washproof, compact and reusable
- Macromolecule foam keeps splint clean and dry



98232	Splint - Finger - 5 x 10cm	B
98234	Splint - Leg - 11 x 100cm	B
98236	Splint - Arm - 11 x 50cm	B

Eye Wash Stations and Solution

Eye Pod Station - Wall Panels

Wall mountable Eye Pod Station with easy to use pods.



80260	Eye Pod Station Mirror, 10 x 20ml Eye Pods & 2 x No.16 Eye Pads	A
-------	--	---

Eye Wash Solution Refills Sterile Saline Solution



Sterile Saline Solution in 500ml bottles, cap designed for quick and easy access to flush the eye thoroughly.

Sterile Saline Solution in 20ml pods, perfect for a smaller incident where a 500ml bottle isn't required.

80202	Eye Wash Solution - 20ml - Box of 25	A
80251	Eye Wash Solution - 250ml	A
80204	Eye Wash Solution - 500ml	A

The Health & Safety Regulations State:

"Where tap water is not available, sterile water or sterile normal saline (Sodium Chloride 0.9%), in disposable containers each holding at least 500ml should be kept easily accessible, and near to the First Aid Box, for eye irrigation. At least 1 litre (2 x 500ml) should be provided."

Double Eye Wash Station

Eye Wash Station can be wall mounted, quick access ideal for use in emergency.



80230	Eye Wash Station Mirror, 2 x 500ml Solution & 2 x No.16 Eye Pads	A
-------	---	---

Combination Wall Panels

The Combination Eye Wash and Pod Station is ideal for smaller incidents where use of 500ml solution isn't necessary.



80248	Combination Eye Wash Station Mirror, 2 x 500ml Solution, 5 x 20ml Pods & 2 x No.16 Eye Pads	A
-------	--	---

Contents	80230	80248	80260
Container	1	1	1
Eye Wash Solution - 500ml	2	2	-
Eye Wash Solution - 20ml Pod	-	5	10
No.16 Sterile Eye Pad	2	2	2
Mirror	1	1	1

Deluxe Eye Wash Cabinet

Deluxe Eye Wash Station complete with bracket which allows for fixing at a specific location. The box can be unclipped from the bracket and taken to the patient.



80226 **Deluxe Eye Wash Station**
Clip in Wall Bracket, 2 x 500ml
Solution and 2 x No.16 Eye Pads

A

Triple Eye Wash Station

Triple Eye Wash Station complete with hanging bracket.



80219 **Triple Eye Wash Station**
Hanging Bracket
3 x 500ml Solution & 2 x No.16 Eye Pads

A

Metal Eye Wash Bracket

The robust metal bracket is suitable for mounting in a vehicle, to hold 1 x 500ml eyewash solution.



80222 **Metal Eye Wash Bracket**
c/w 1 x 500ml Solution

A

80218 **Metal Eye Wash Bracket**
Empty

A

Premium Eye Wash Station

The Premium Eye Wash Station complete with hanging bracket.



34226 **Premium Eye Wash Station**
Hanging Bracket, 2 x 500ml Solution,
5 x 20ml Pods and 2 x No.16 Eye Pads

A

COSHH Regulations State:

"Employers are to provide emergency eye/skin wash facilities wherever the work force may be exposed to injury from chemicals or any substance hazardous to health".

Contents

	80222	80226	80219	34226
Container	1	1	1	1
Eye Wash Solution - 500ml	1	2	3	2
No.16 Sterile Eye Pad	-	2	2	2
Wall Hanging Bracket	-	-	1	1
Clip in Wall Mounting Bracket	-	1	-	-
Mirror	-	-	-	1

Bio Hazard Clean up



77918	Super Absorbent Granules - 100gm	C
77922	Super Absorbent Granules - 500gm	C
77932	Clean Up Pack With Oversleeves	C
77906	Disinfectant Spray 500ml (clear)	C

First Aider Protection Packs

These packs help provide protection against cross infection for first aiders. Ideal for isolated staff members e.g. school playground staff, security personnel and warehouse operatives.

91810	Pack A - Personal Protection Min-Grip Poly Bag	B
P91810	Pack B - PU Nylon Pouch With Belt Loop c/w 1 x Content	B
P91812	Pack C - PU Nylon Pouch With Belt Loop c/w 2 x Content	B

Bio Hazard Combination Kit

- 2 x Body Fluid Spill Kit Refills
- 1 x Sharps Disposal Kit c/w Disinfectant
- Supplied in Premium Bio Hazard Box



70720	Combination Kit	B
-------	-----------------	---

Body Fluid Spill Kit

Body Fluid Spill Kits – 1 application content as follows:

- Disinfectant: a broad spectrum viricide, bactericide and fungicide, highly effective against micro-organisms (Viruses, bacteria and fungi)
- Super absorbent: a fast acting, sugar absorbent powder that rapidly congeals and deodorises fluid spillages
- Selected clothing: polythene oversleeves, aprons, disposable gloves, bio hazard bag, plastic scoop and scrape
- Accessories: non-woven wipes and an antiseptic wipe

77904	Body Fluid Refill Pack 1 application	B
70712	Body Fluid Two Application Kit Size 227 x 165 x 47mm	B
70714	Body Fluid Five Application Kit Size 227 x 165 x 47mm	B

Contents: (single use only items)

- 1 x Resusciate Face Shield
- 1 x Disposable PVC Apron
- 1 x Disposable Vinyl Gloves
- 1 x Cleansing Wipe
- 1 x Bio Hazard Disposal Bag

IDEAL
TO INCLUDE
IN A FIRST
AID KIT



Wipes



82124	Moist Cleansing Wipes - Box 100	B
82122	Isopropyl Alcohol Wipes - Box 100	B

Sharps Handling & Disposal

This range is designed for the potentially hazardous task of handling discarded sharps such as needles and syringes, blades and other sharp objects that may be contaminated.

IMPORTANT:

- Always assess the elements of risk before commencing the removal task
- Always make sure that all protective clothing and equipment required is available
- Never put any part of the body into areas where there is a restricted or unclear view

Hex Armor Gloves

Highest needlestick resistance available provided by the layering of SuperFabric brand material

- Tested using actual 25 gauge hypodermic needles
- Highest level of cut resistance in enhanced areas; European standard EN388 Cut Level 5

Sharps Disposal Kits

Sharps such as contaminated needles and syringes, knives, blades and broken glass can puncture the skin and act as a vehicle for the transmission of harmful micro-organisms.



96020	Sharps Disposal Kit c/w Disinfectant Spray	B
96021	Sharps Disposal Kit	B

Sharps Handling Grabber



96080	Sharp Safe Mechanical Grabber	C
-------	-------------------------------	---

Sharps Handling Kits



96030	Sharps Handling Kit c/w 1 Pair Of Hex-Armour Glove	D
96032	Sharps Handling Kit Without Glove	D
96090	Hex Armour Needle Protection Glove 1 pair	D

Sharps Disposal Containers



96000	Sharps Container - 0.6 Litre	B
96002	Sharps Container - 2 Litre	B
96006	Sharps Container - 4 Litre	B
96010	Sharps Container - 7 Litre	B
96012	Sharps Container - 24 litre	B
96014	Sharps Container - 30 litre	B

First Aid Room Packages

Economy First Aid Room Furniture Package

- Examination couch with paper roll holder
- Single tier step
- 3 Panel mobile folding screen
- Large stainless steel first aid cabinet (1 shelf)
- 46cm wide dressing trolley with two glass effect safety trays
- 20 Litre Bin
- Delivery normally within 3-4 weeks



Standard First Aid Room Furniture Package

- Deluxe Examination couch with paper roll holder
- Single couch drawer unit in white
- Single tier step
- 4 Panel mobile folding screen
- 46cm wide dressing trolley with three glass effect safety trays and drawer unit in white
- Medium stainless steel first aid cabinet (2 shelves)
- 20 Litre Bin
- Delivery normally within 3-4 weeks



Premium First Aid Room Furniture Package

- Deluxe Examination couch with paper roll holder
- Three couch drawer units in white
- Single tier step
- 4 Panel mobile folding screen
- Trolley with 2 single and 2 double clear trays
- 50cm Wide locking base cabinet - white
- 50cm Wide locking wall cabinet - white
- 20 Litre Bin - white
- Delivery normally within 3-4 weeks



Deluxe First Aid Room Furniture Package

- Specialist treatment couch with paper roll holder
- One drawer line unit & one drawer pack in white
- Double tier step
- 4 Panel mobile folding screen
- Trolley with 4 single and 1 double clear tray
- 100cm Wide locking base cabinet - white
- 100cm Wide locking wall cabinet - white
- 50 Litre Bin - white
- Delivery normally within 3-4 weeks



School First Aid Room Furniture Package

- Rest couch with paper roll holder
- 4 Panel mobile folding screen
- Trolley with 4 single blue trays
- 50cm Wide locking base cabinet - white
- 50cm Wide locking wall cabinet - white
- 20 Litre Bin - white
- Delivery normally within 3-4 weeks



First Aid Room

Column Scales

Features:

- Eye-level reading
- Low profile footplate
- Built-in height measurement



95064 Column Scales - 149 x 43 x 53cm

B

Pedal Bin

Safely and hygienically store clinical waste ready for disposal.

Features:

- Removable inner plastic bucket with handle
- 20 litre capacity
- Easy clean surface and interior

Specification:

- Size: 45cm (H) x 30.5cm (D)
- Weight: 2.8kg



96854 Pedal Bin - 20 litre

B

Trolley

HSE recommended for all First Aid rooms.

Features:

- Three 'lipped' shelves for larger item storage.
- Castored feet with locking pedals.
- Easy to clean

Specification:

- Size: 90cm (H) x 62cm (W) x 45cm (D)
- Weight: 12.5kg



95160 Medical Trolley - 3 Shelves

B

Step

Lightweight step for treatment couch.

Features:

- Tubular steel frame
- Non-slip feet
- Easy to clean

Specification:

- Size: 30cm (H) x 40cm (W) x 30cm (D)
- Step Weight: 2kg
- Maximum Load: up to 159kg



95111 Step

B

Medical Screen

With four-fold curtain design.

Features:

- Stainless steel frame
- Removable curtains
- Simple manoeuvrability
- Folds away for neater storage

Specification:

- Size: Open 179.5cm (H) x 235cm (W)
- Size: Closed 179.5cm (H) x 62.5cm (W)
- Weight: 12.5kg



95316 Four Panel Medical Screen - Blue

B

Portable Couch

The couch transforms from a convenient sized carry case weighing less than 20kg into a full sized couch.

Features:

- Adjustable heights
- Padded armrests
- Cut out face panel
- Complete with carry case
- Black and powder coated metal frame
- Easy to clean

Specification:

- Size with attachments: Open 88cm (H) x 96cm (W) x 216cm (L)
- Size without attachments: Open 88cm (H) x 66cm (W) x 188cm (L)
- Size: Folded 16cm (H) x 70cm (W) x 95cm (L)
- Couch Weight: 18.5kg
- Maximum Load: up to 159kg



Treatment Couch with Adjustable Backrest

An ideal item in the workplace, GP Surgery and most medical environments the couch is simple to maintain, easy to clean and built to last.

Features:

- Couch roll holder
- Fixed height
- Two sections
- Adjustable head/back rest from horizontal to 85 degrees
- Soft faux leather upholstery
- Easy to clean and maintain

Specification:

- Size: 78cm (H) x 60cm (W) x 193cm (L)
- Couch Weight: 20.5kg
- Maximum Load: up to 150kg



95110	Treatment Couch w/ Adjustable Backrest	C
95050	Paper Roll - 20m	C

EVACUATION & RESCUE

Deluxe Evacuation Chair

Features:

- Special track structure to aid with ascending and descending stairs
- Space saving, can be folded completely during storage
- Multiple restraint straps with quick release buckles
- High strength aluminium material
- Four wheels & Footrest

Specification:

- Size: 70 (L) x 54 (W) x 92cm (H)
- Folded size: 70 (L) x 54 (W) x 17cm (H)
- Load bearing: 159kg
- Weight: 8kg



95081

Deluxe Evacuation Chair

C

Mobile Transit Chair

Features:

- Adjustable handles aiding patient transportation
- Light weight & Small compact size
- Easy to use and manoeuvre
- Easy to clean and sterilise

Specification:

- Size: 90 (L) x 51 (W) x 91cm (H)
- Folded size: 90 (L) x 17 (W) x 59cm (H)
- Load bearing: 159kg
- Weight: 8kg
- Colour: Blue



95082

Mobile Transit Chair

C

Evacuation Chair

Our Evacuation Chair features an advanced tracked design allowing use by one operator, reducing health and safety concerns during use and allowing smooth decent of stairs.



Features:

- Complete with cover and wall bracket
- Lightweight and compact for easy storage
- Adjustable head pad and patient restraint
- Impervious to bodily fluids
- Aluminium alloy materials
- Easy to clean
- 5 year warranty

Specification:

- Size Open: 136cm (H) x 52cm (W) x 90cm (D)
- Size Folded: 112cm (H) x 52cm (W) x 20cm (D)
- Chair Weight: 10kg
- Maximum Load: up to 159kg

95080

Evacuation Chair

C

Patient Wheelchair

Ideal as a part time wheelchair, these chairs also prove very popular in nursing homes, hospitals and public areas.



Features:

- Folding back
- Fixed arm rests
- PU wheelchair tyres
- Storage pouch on back rest
- Rear wheel locking brakes
- Neatly folds for storage or transportation
- Swing away foot rests with height adjustment

Specification:

- Maximum User Weight: 115kg or 18 Stone
- Device Weight: 14.5kg or 32lbs
- Overall Height: 89cm • Height Folded: 74cm
- Overall Width: 59.5cm • Width Folded: 30cm
- Overall length: 93cm • Length without Footrests: 46cm

95025 Wheelchair

C

Scoop Stretcher

Firm, lightweight, easy to clean and operate.



Features:

- Twin safety lock for smooth locking and unlocking
- Recessed head section to maintain patient's cervical alignment
- Four length options to accommodate the patients height
- Thermally-treated polymer material
- X-ray translucent

Specification:

- Product Size: 7cm (H) x 44.5cm (W) x 162cm (L)
- Gross Weight: 8kg • Loading Bearing: 159kg

98204 Break Away Scoop Stretcher

C

Neil Robertson Stretcher

Features:

- Designed for confined space rescues • Folds neatly during storage
- Multiple patient restraint straps keep patient immobilised
- Ideal for situations requiring vertical lifting
- Used by the Royal Navy, Merchant Navy, in Mining Industries and the Emergency Services.

Specification:

- Size: 165 (L) x 102 (W) x 5cm (H)
- Load bearing: 159kg & Weight: 3.5kg



98202

Neil Robertson Stretcher

C

Vacuum Immobilisation Mattress

Features:

- Capable of immobilising patients with all kinds of injuries including spinal and neck
- Constructed from isolative and cold resistant PVC
- Compatible with X-ray and MRI examinations
- Promotes patient comfort and reduces the risk of pressure sores.
- Safely adjust the stretcher to be soft or hard using the dual action pump.

Specification:

- Head width: 80cm
- Middle width: 88cm
- Feet width: 80cm
- Length: 200cm
- Weight 10kg



98102

Vacuum Immobilisation Mattress

C

Evacuation & Rescue

Duo-Fold Stretcher

The Duo-Fold stretcher requires minimum storage space when not in use, folds both horizontally and vertically into the compact dimension shown. The frame is manufactured from aluminium and is both strong and lightweight.



Features:

- High strength aluminium alloy frame
- Waterproof, fireproof and anti-cracking
- Complete with fully zipped carry bag and straps
- Hi-Vis Orange
- Easy to clean

Specification:

- Size: Open 14cm (H) 52cm (W) 206cm (L)
- Size: Folded 10cm (H) x 18cm (W) x 103cm (L)
- Stretcher Weight: 5.5kg
- Maximum Load: up to 159kg

98106

Duo Fold Pole Stretcher

C

Basket Stretcher

Sturdy, flexible and is used in all types of rescue/recovery situations. With its special lifting sling equipment, the stretcher is ideal for use with a helicopter winch.



Features:

- Adjustable feet, safety belts and mattress
- Fire resistant
- High resistance to everyday wear and corrosion

Specification:

- Product Size: 19cm (H) x 61cm (W) 216cm (L)
- Load bearing: less than 159kg
- Gross Weight: 20kg

98227

Basket Stretcher

C

98228

Basket Stretcher Lifting Harness

C

Basket Stretcher - 2 Piece

Sturdy, flexible and is used in all types of rescue/recovery situations. With its special lifting sling equipment, the stretcher is ideal for use with a helicopter winch.



Features:

- Easy to store, breaks apart into 2 pieces
- Adjustable feet, safety belts and mattress
- High resistance to everyday wear and corrosion

Specification:

- Product Size: 216 (L) x 61 (W) x 19cm (H)
- Folded Size: 130 (L) x 61 (W) x 26cm (H)
- Packing size: 135 x 64 x 28cm
- Load bearing: less than 159kg
- Weight: 20.5kg
- Fire resistant

98225

Basket Stretcher - 2 Piece

C

Spinal Board

The spinal board is made of PE material with no discharge contaminator, and firm to wear. Typically used to carry patients and injured persons after sustaining injury outdoors.



Features:

- Floatable device for water rescue situations
- Compatible with most head immobilisation devices and strap mechanisms
- X-ray translucent

Specification:

- Product size: 7cm (H) x 45cm (W) x 184cm (L)
- Load bearing: less than 159kg
- Gross Weight: 9kg

98206	Spinal Board	C
-------	--------------	---

Patient Transfer Sheet



Features:

- Lightweight and compact, supplied with carry case
- Multiple carry straps with ten strategically placed handles
- Easy to clean vinyl material

Specification:

- Size: Open 204cm (L) x 71cm (W)
- Maximum Load: up to 159kg

98229	Patient Transfer Sheet c/w Carry Case	C
-------	---------------------------------------	---

Cervical Collar



The foam cervical collar is used for cervical injuries requiring non-rigid support and limited mobility. Can also be used as an aid for the relief of pain associated with neck muscle or whiplash injuries.


Features:

- Soft foam is slightly contoured for a comfortable fit
- Covered with 100% cotton stockinet and a Velcro closure
- Loop contact closures for adjustability

98240	Cervical Collar - Blue	C
-------	------------------------	---

98241	Cervical Collar - Grey	C
-------	------------------------	---

Head Immobiliser



Made of high-density plastic material, large ear holes for monitoring the patient's ear canal.

Features:

- Compatible with scoop stretcher and spine board
- Minimal interference with X-Ray, MRI, or CT scanning procedures
- Water-resistant plastic coating makes it easy to clean and prevents bacterial growth

98217	Head Immobiliser	C
-------	------------------	---

Defibrillators



Mediana Semi Automatic Defib

- Supplied with a protective outer cover
- Switches on automatically when opened
- Defibrillator pads are constantly connected and ready for use
- A combination of LED indicator lights, visual icons and clearly spoken audible voice prompts, make the AED safe to use even by untrained persons
- Simple flip switch sets the machine to either Adult or Paediatric, the pads are suitable for both situations
- The battery is a disposable LiMnO2 Non-Rechargeable version with two years shelf life, and five years standby life

95604

Mediana Semi-Automatic Defibrillator

B

Mediana Semi Automatic Defib Bundle

Includes:

- Mediana Semi Automatic Defib
- SD Card
- AED Prep Kit

AED Prep Kit contents:

- 1 Cleansing Wipe
- 1 Clinical Waste Bag
- 1 Disposable Razor
- 1 Gloves (Pair)
- 2 Hand Towels
- 1 Resusciate
- 1 Heavy Duty Clothing Shears



95624

Mediana Defib Bundle Pack

B

Pads

- Pads are suitable for both Adult and Paediatric



95626

Mediana Defib Replacement Pad Pack

B

Battery Pack

- Disposable
- LiMnO2 Non-Rechargeable
- 2yrs shelf life
- 5yrs standby life



95625

Mediana Defib Battery Pack

B



Outdoor AED Cabinet - Alarmed

Specification:

- Power 200W/230W
- Adjustable Bulb Thermostat -30°C to 30°C
- Armour Resistance AM
- d.8/22, AISI321, 1x200W, 23V
- Cable 1.3m
- Degree of Protection IP44
- Installation requires hard wiring to suitable outlet

Features:

- Manufactured from 4mm thick Carbon Steel
- Shatterproof Transparent Window
- Outdoor Weather Resistance
- Anti-theft Radio with Magnetic Sensor
- Thermally Regulated Heating System
- Alarm System (105db) with battery back up
- Independent remote to arm/disarm the siren
- Magnetic sensor
- Optional padlock
- Safe and Visible
- Wall Mountable

95606	AED Outdoor Heated Wall Cabinet Alarmed - 480 x 425 x 215mm	B
95605	AED Metal Wall Cabinet Alarmed - 460 x 300 x 140mm	B

AED Metal Wall Cabinet - Alarmed

- Toughened durable glass door
- Compatible with most AED models
- Alarm sounds when opened
- Lockable
- Safe and Visible
- Quick Release Push-Button
- Fixed shelf for added storage
- Wall mountable
- Size 46cm (H) x 30cm (W) x 14cm (D)



Universal AED/Defibrillator Bracket


- Sturdy metal bracket is ideal for mounting any AED on a wall
- 12.3cm (H) x 19cm (W) x 9.5cm (D)



95607	Universal AED/Defib Wall Bracket	B
-------	---	---

Universal AED/Defibrillator Sign

- Self-Adhesive Vinyl
- 200 x 150mm



6053/S-A	Defibrillator Sign - 200x150 Self Adhesive	B
----------	---	---

First Aid Equipment and Accessories

Resuscitators

- Latex free medical grade PVC with polycarbonate shell
- Clear mask and valve housing
- Anatomically shaped mask with soft adjustable air cushions
- Superb seal and patient comfort
- Infant and Child sets feature a 40cm H2O pressure relief valve
- Supplied bagged and ready for use
- Single use to help to prevent cross contamination



Resuscitation



91808	Resuscitade Face Shield	B
91815	Resuscitade Face Shield c/w Keyring Pouch	B
91805	Pocket CPR Face Mask c/w Plastic Case	B

Guedel Airways

- Rigid bite block to ensure airway patency
- Smooth integrated design for patient comfort and safety
- Colour coded for easy selection of size
- Sterile individually wrapped
- Full range of sizes available



99710	Guedel Airway *Prefix code with size required	B
-------	--	---

Thermometers



91130	Office & Shop Thermometer	C
91132	Clinical Ear Thermometer Covers	C
91134	Clinical Ear Thermometer	C
91135	Disposable Thermometers - 100 pk	C
91136	Digital Thermometer	C

Thermometers

95002	Digital Blood Pressure Monitor	C
-------	--------------------------------	---

- Compact and easy to use
- Accurate semiconductor sensor with low noise air pump
- Integrated wrist cuff



Tweezers



91118	Tweezers - 8cm Nickel Plated	B
91122	Tweezers - 15cm Blunt End	B
91124	Tweezers - 13cm Fine Point	B
91126	Tweezers - Plastic	B

Trauma Foil Blanket

Orve+wrap® creates a convection heat to ensure a hypothermic or trauma patient is pro-actively warmed and will maintain normothermia

91008	Orve+Wrap 200 x 150cm	B
-------	--------------------------	---

Emergency Foil Blankets



91004	Emergency Foil Blanket - 210 x 160cm	A
-------	--------------------------------------	---

Syringes

99016	Terumo 1ml 3 Part L/S Syringe - Box 100	-
99017	Terumo 2ml 3 Part L/S Syringe - Box 100	-
99021	Terumo 5ml 3 Part L/S Syringe - Box 100	-
99023	Terumo 10ml 3 Part L/S Syringe - Box 100	-
99030	Terumo 20ml 3 Part L/S Syringe - Box 50	-
99029	Terumo 30ml 3 Part L/S Syringe - Box 50	-

Scissors



91100	Scissors - Sharp / Sharp	B
91102	Scissors - Blunt / Blunt	B
91104	Scissors - Blunt / Sharp	B
91112	Heavy Duty First Aid Snips	B

First Aid Blanket



91022	First Aid Blanket - Cotton Cellular 150 x 200cm	B
-------	--	---

Misc. First Aid Equipment

91424	Accident Record Book	B
90098	First Aid Armband Elasticated	B
91502	First Aid Enamel Badge	B
97024	Pen Torch	B
96840	Kidney Dish	B
95000	Stethoscope	B
80254	Plastic Eye Bath	B
97026	Wooden Tongue Depressors - Box 100	B
71196	Security Seals - 10 pk	B

Fire Extinguishers

Fire Risk Assessment Guide

Commissioning of a new extinguisher and servicing should be carried out by a trained person as set out in British Standard BS: 5306-3: 2009.

- ✓ **Risk Assessment**
Is there an appropriate recording system?
- ✓ **Fire Fighting and Fire Detection**
Is there appropriate equipment on site?
Is fire fighting equipment clearly identified and locatable?
Are there nominated people who are trained for fire emergencies?
Are all emergency routes and exits clearly identified and illuminated?
- ✓ **Emergency Escape Routes**
Do emergency routes or exits lead to a place of safety?
Are everyone's needs catered for in event of an emergency escape?
Are dangerous areas restricted to authorised personnel only?
- ✓ **Maintenance**
Is fire fighting equipment fully maintained and serviced with written documentation?

Under The Regulatory Reform (Fire Safety) Order 2005, the employer or 'responsible person' for that building must ensure that a fire risk assessment has been conducted and that appropriate steps have been put into place to remove or reduce risks.

- ✓ **Procedures for Serious/Imminent Dangers and Areas**
Are all authorised persons given equipment to protect them from risks?
Are all dangerous substances clearly identified and safety data sheets provided?
Are all dangerous substances kept in appropriate storage?
- ✓ **Provision of Information for Employees**
Can the responsible person provide comprehensive information of risks to employees as identified by the risk assessment?
Are all off-site employees given relevant information on fire risks and preventive and protective measures on the premises?
- ✓ **Training**
When they began their employment, were all employees given adequate safety training?
Has all training been officially recorded?

Fire Classification	Foam	ABC Powder	CO2	Water	Wet Chemical
Flammable solids Wood, paper, textiles, fabrics etc.	YES Smothering & cooling effect	YES Rapid smothering effect prevents re-ignition	NO	YES Saturation & cooling prevents re-ignition	YES
Flammable liquids Petrol, oil, fat, paint, solvents, grease etc.	YES Blanket fires preventing re-ignition	YES Rapid results	YES Rapid results leaves no residue	NO	YES Fat Fires
Flammable gases Propane, Butane, Natural gases	NO	YES	NO	NO	NO
Electrical equipment Computers, printers, etc.	NO Some restricted use	YES	YES Non-conductive	NO Water is conductive	YES 35KV dielectric test

Powder Fire Extinguishers

MediKit's range of Dry Powder extinguishers are suitable for class A, B and C fires, meaning they can tackle most fires involving wood/paper, flammable liquids and electrical fires.

WG10106	600GM ABC 5A/21B	C
WG10102	800GM ABC 21B	C
WG10114	1.0KG ABC 8A/34B	C
WG10124	2.0KG ABC 13A/89B	C
WG10132	3.0KG ABC 21A/113B	C
WG10134	4.0KG ABC 27A/144B	C



WG10142	6.0KG ABC 43A/223B	C
WG10148	9.0KG ABC 55A/233B	C



Foam Fire Extinguishers

Our foam extinguishers achieve high fire ratings for both A and B class fires, making them highly effective against ordinary wood/paper fires (A class) and also petrol and volatile liquids (B class). The foam forms a flame smothering seal over the surface of the burning material which prevents re-ignition.

- 35kV dielectric test
- Internal Polyethylene lining
- Complete with spray nozzle

WG10112	1.0LTR 5A/21B	C
WG10130	2.0LTR 5A/34B	C
WG10144	6.0LTR 21A/144B	C
WG10149	9.0LTR 27A/233B	C

Wet Chemical Fire Extinguishers

Designed for Class F fires, Wet Chemical extinguishers are ideal for kitchens and caravans where cooking oils and fats are the risk. Independent fire tests have established that it is capable of extinguishing a fire in a 25 litre deep fat fryer.



WG10146	6.0LTR 21A/113B/75F	C
---------	---------------------	---

CO2 Fire Extinguishers

CO2 extinguishers are used for situations involving flammable liquids and electrical hazards. The CO2 extinguishant is harmless to delicate equipment making it ideal for office environments.

- Harmless to machinery
- Aluminium alloy construction



WG10128	2.0KG 34B	C
WG10138	5.0KG 70B	C

Water Fire Extinguishers

Designed, developed and tested in the UK for use on A-class fires, this product carries the well-known Kitemark, LPCB and MED certification achieving superior results in tests.

- Includes 35kV dielectric test
- Fluoro Free additive
- Internal Polyethylene lining



WG10140	6.0KG 13A	C
WG10150	9.0KG 21A	C

Vehicle and Fire Safety

BEST SELLER



Vehicle Fire Extinguishers

This range is effective on vehicle fires making them great extinguishers to have in the car or van. The premium travel brackets have a clip in metal arm which hold the extinguisher securely in place during transit.

- Fully Approved by LPCB (UKAS Accredited)
- Approved to BSEN3
- Holds 5 year guarantee

PREMIUM TRAVEL BRACKET

WG10116	1.0KG ABC 5A/34B Premium Travel Bracket	B
WG10126	2.0KG ABC 8A/55B Premium Travel Bracket	B

Vehicle Fire Extinguishers

The small size of this extinguisher makes it ideal for cars, vans, caravans, taxis and the home.

- EN3 Approved (except clause 6.1 size)
- UKAS accredited
- Holds 5 year guarantee



Advantages of dry powder

- Quick knockdown when correctly applied
- Safe to use on Electricity up to 1000 volts
- Effective on running fires such as petrol, diesel and oils
- Creates no thermal shock for fires onto hot metal, effective on undercarriage fires
- Good heat shield, creates an effective shield against radiated heat
- Can be used on all class fires A, B and C
- Effective on a range of temperatures from 15 to +60 degrees Celsius, unlike water based extinguishers such as foam, which can freeze.

Vehicle Safety Kit

Suitable for Private Motorists supplied in a zip bag complete with handle/shoulder strap.

Contents:

- One Person Travelling First Aid Kit
- Fire Extinguisher
- Warning Triangle
- Glass Hammer



Winter Safety Kit

An essential kit for the winter months, supplied in a nylon bag for easy storage within a vehicle.

Contents:

- Folding Snow Shovel in Nylon Pouch
- One Person Travelling First Aid Kit
- Foil Blanket
- Hi-Vis Vest



70147	Vehicle First Aid Kit Soft Nylon Bag - 470 x 150 x 140mm	A
--------------	--	----------

33142	Winter Motoring Kit Soft Nylon Bag - 470 x 150 x 140mm	A
--------------	--	----------

DIN Euro Motoring Kit

Complies with DIN 13164 European Specification.

Contents:

- 1 Burn Dressing (Sterile) 40 x 60cm
- 1 Burn Dressing (Sterile) 60cm x 80cm
- 1 Standard Dressing (Sterile) 10cm x 12cm
- 1 Standard Dressing (Sterile) 6cm x 8cm
- 2 Standard Dressings (Sterile) 8cm x 10cm
- 2 Triangular Bandages 136cm x 96cm x 96cm
- 3 Wound Compresses 10cm x 10cm (2pcs)
- 14 Plasters, Washproof Assorted
- 2 Conforming Bandage 6cm x 4m
- 3 Conforming Bandage 8cm x 4m
- 1 Adhesive Tape 2.5cm x 5m
- 2 Cleansing Wipes
- 2 Gloves (Pairs)
- 1 First Aid Leaflet
- 1 Scissors



48464	Hi-viz Yellow Waistcoat/Vest Sizes - XXL, XXXL	B
41996	Travel Sickness Disposal Bag - Pack 12	B
66020	Vehicle Warning Triangle	B
66026	Vehicle Seat Belt Cutter	B
66028	Vehicle Safety Hammer for Glass	B
66050	Folding Shovel	B
97020	Snap-Light - Single Use Stick	B
97084	NF Alcohol Breath Test - Twin Pack	B
91004	Emergency Foil Blanket	A

70144	DIN 13164 Euro Motoring Kit	A
-------	-----------------------------	---

Fire Marshal Kits

These comprehensive kits contain Emergency Fire equipment suitable for any workplace. Keep all Emergency Fire equipment in one practical bag, there's even room for extra supplies.

Standard Fire Marshal Kit

Contents:

- Fire Warden Vest
- Reflective Armband
- 1 x Airhorn
- 6 x Plastic Whistles
- 1 x Torch
- 10 x Light Sticks
- 1 x Mini Megaphone



Deluxe Fire Marshal Kit

Contents:

- 1 x Small BS 8599-1 First Aid Kit
- 1 x Small Burns Kit
- 6 x Emergency Foil Blankets
- Fire Warden Vest
- Reflective Armband
- 1 x Airhorn
- 6 x Plastic Whistles
- 1 x Torch
- 10 x Light Sticks
- 1 x Mini Megaphone
- 6 x Nitrile Gloves



WG30500	Fire Marshal Kit 410 x 350 x 285mm	A
---------	---------------------------------------	---

WG30502	Deluxe Fire Marshal Kit 435 x 420 x 360mm	A
---------	--	---

Fire Safety Products

Fire Warden Vests

- Fire Warden Emblem transferred
- EN471 Class 2 compliant
- Three sizes available



WG30106	Fire Warden Vest XL	C
WG30108	Fire Warden Vest M	C
WG30110	Fire Warden Vest L	C

Pack-a-Cone

- Pull-up traffic cone that can be deployed and packed away in seconds
- Ideal for emergencies and breakdowns



WG30212	Pack a Cone	C
---------	-------------	---

Megaphones

- Incorporates integrated siren for emergency situations
- Rugged ABS plastic construction
- Press to talk trigger



WG30300	Mini Megaphone	C
WG30302	Megaphone	C

Fire Blankets

- Kitemark certified to BSEN1869:1997
- Front loading design for minimal removal time
- Tamper evident security feature
- Requires 40% less wall space



WG30400	IVG Fire Blanket - Plastic Case	C
WG30402	Fire Blanket - Soft Pouch	C

Smoke Alarms


- Fully certified to EN14604
- Supplied complete with battery
- Test button verifies battery and alarm operation
- Low battery indicator



WG30206	Smoke Alarm w/ Light	C
WG30208	Micro Smoke Alarm	C

Keybox with hammer

- Easy to use 'snatch off' design
- Spare lenses available



WG10234	Plastic Keybox and Hammer	C
---------	---------------------------	---

Fire Log Book

- Central record for all aspects of fire safety procedures
- Includes testing and maintenance guidelines for each section



WG30314	Fire Log Book - Pack 10	C
---------	-------------------------	---

Box

Top-loading cabinets with internal spring mechanism designed to hold the extinguisher in place whilst in transport.

- 6kg/9kg Cabinet
- Capacity: 1 x 6kg/9kg extinguisher
- Totally corrosion-proof design
- Neoprene seal
- Polypropylene vision panels in door
- Internal equipment retaining straps
- Available with integral alarm



WG10230	6.0KG/9.0KG Vehicle Cabinet	C
WG10224	6.0KG Vehicle Box	C
WG10225	9.0KG Vehicle Box	C

Air horn

- Loud, long sound duration (120dB)
- Single, non-flammable gas propellant
- Ideal for providing temporary cover during alarm installations



WG30304	Airhorn 300ml	C
---------	---------------	---

Stands

- Holds 1 extinguisher up to 12kg
- Non-GRP construction
- Suitable for use in food prep areas
- Features skirting board cut out



WG10220	Fire Extinguisher Stand Single	C
WG10222	Fire Extinguisher Stand Double	C

First Aid Boxes & Bags

ABS Empty First Aid Boxes

MediKit's Deluxe First Aid box is a tough and durable ABS plastic case complete with wall bracket and internal compartments.

- The box can be opened whilst fixed in the bracket without spilling the content
- These boxes also have rubber gasket seals which help prevent contamination of contents by dust and moisture



AB3	ABS Box - 270 x 185 x 115mm	A
AB4	ABS Box - 320 x 215 x 125mm	A
AB5	ABS Box - 400 x 300 x 150mm	B

Plastic Boxes

- Tough and durable plastic box
- Optional clip in wall



P6	Plastic Box - 240 x 170 x 80mm	A
P8	Plastic Box - 270 x 220 x 85mm	A
P9	Plastic Box - 340 x 270 x 100mm	A
32156	Clip In Wall Bracket - to take P6	B
32157	Clip In Wall Bracket - to take P8	B
32158	Clip In Wall Bracket - to take P9	B

Empty First Aid Bags



D	Medium Green Zip Pouch c/w Belt Loop 180 x 110 x 45mm	A
D1	Large Green Zip Pouch c/w Belt Loop 210 x 140 x 60mm	A
B1	Medium Velcro Belt Pouch c/w Belt Loop 170 x 190mm	A

- Durable, lightweight fabric
- Suitable for heavy duty clothing shears (Tufkut) or a torch
- Complete with belt loop

B12	Accessory Belt Pouch 200 x 50 x 40mm	A
-----	---	---

B3	Large Green Velcro Belt Pouch c/w Belt Loop 190 x 165mm	B
D5	Deluxe Green Zip First Aid Pouch c/w Belt Loop 255 x 155 x 95mm	B
C2	First Aid Green Bum Bag c/w Waist Strap 420 x 125 x 120mm	B



First Aid Holdall

- Multi compartments
- Detachable shoulder straps
- Internal compartments for easy storage and access
- The bags have extra space which can be used to include additional products



First Aid Grab Bag

- Detachable shoulder straps
- Internal compartments for easy storage and access
- The bags have extra space which can be used to include additional products



SB7 **Large Green First Aid Carry Bag**
360 x 280 x 180mm
Shoulder Strap and Internal Compartments

A

A **Green First Aid Shoulder Bag**
280 x 230 x 120mm
Shoulder Strap and Internal Compartments

A

Active Bag Range

Our Active Bag range has been designed with outdoor pursuits in mind: durable, rip-resistant and easy to clean.



The compartments make the Active Holdall Bag range practical in all environments with easy organisation of supplies, even providing room for extra content when bought complete.

- Combination of up to 14 Zip and Velcro compartments plus additional storage
- Protective padded lining
- Reflective strips
- Strong, water-resistant, highly visible material

RS3 **Holdall - Green - 260 x 380 x 170mm**

A

RS4 **Holdall - Orange - 260 x 380 x 170mm**

A



Super compact and versatile, can be used as a shoulder bag, bum bag or simply carried by hand. Durable, rip resistant and easy to clean our Active Bag is extremely practical in all pursuits.

- Internal and adjustable divider
- Internal, external and side pockets
- Adjustable shoulder and belt straps
- Carry handle
- Safety clasps
- Secure key holder
- Strong, water-resistant, highly visible material

RS1 **Bum Bag - Green- 175 x 300 x 100mm**

A

RS2 **Bum Bag - Orange - 175 x 300 x 100mm**

A

First Aid Boxes & Bags

Trauma Barrel Bag

- Trauma, fold out multi compartment barrel bag
- Durable, water and abrasion resistant material
- Padded sections to carry oxygen cylinders
- Versatile bag supplied with carry handle and adjustable shoulder strap
- Reflective strips



PR72 **Trauma Barrel Bag**
495 x 170 x 260mm

A

Medical Rucksack

- Fold out, multi compartment trauma rucksack
- Durable, water and abrasion resistant material
- Supplied with carry handle and padded adjustable Rucksack straps with adjusters
- Padded back panel for user comfort



PR8 **Rucksack**
380 x 315 x 145mm

A

Multi Pouch Bag

- Durable, water proof, light weight material
- Abrasion resistant material
- Padded carry handle
- Five inner velcro removable zip clear front pouches
- Super practical bag. contents can be organised specific to use



PR9 **Emergency Bag**
400 x 225 x 95mm

A

Medi Pouch

- Fold out zip pouch with belt loop and carry handle
- Durable, water and abrasion resistant material
- Four compartments, practical organisation of supplies



PR10B **Pouch - Blue - 205 x 120 x 55mm**

A

PR10G **Pouch - Green - 205 x 120 x 55mm**

A

PR10R **Pouch - Red - 205 x 120 x 55mm**

A

Clear Medi Pouch

- Zip medi case with carry handle
- Clear front, easily view medicaments inside
- Durable, water and abrasion resistant material



PR11Y **Clear Medi Pouch - Yellow**
215 x 135 x 90mm

A

PR11G **Clear Medi Pouch - Green**
215 x 135 x 90mm

A

PR11R **Clear Medi Pouch - Red**
215 x 135 x 90mm

A

Extra Large Trauma Bag

- Multi compartment extra large trauma bag
- Various sized compartments for practical organisation of supplies
- Durable, water proof, light weight material
- Abrasion resistant material
- Super versatile supplied with carry handle, adjustable shoulder strap with grip and adjustable back pack straps
- Removable wipe clean base with grip studs
- Reflective strips



PR7	XL Large Trauma Bag 435 x 420 x 360mm	A
-----	---	---

Trauma Barrel Bag

- Trauma, fold out multi compartment barrel bag
- Durable, water proof, light weight material
- Abrasion resistant material
- Padded sections to carry oxygen cylinders
- Versatile bag supplied with carry handle and adjustable shoulder strap
- Reflective strips



PR71	Trauma Bag Wipe Down 495 x 170 x 260mm	A
------	--	---

Large Trauma Bag

- Multi compartment large trauma bag
- Various sized compartments for practical organisation of supplies
- Durable, water and abrasion resistant material
- Super versatile supplied with carry handle, adjustable shoulder strap with grip and adjustable back pack straps
- Removable wipe clean base with grip studs
- Reflective strips



PR6	Large Trauma Bag - Red 410 x 350 x 285mm	A
-----	--	---

PR6G	Large Trauma Bag - Green 410 x 350 x 285mm	A
------	--	---

Fold Out Paramedic Bum Bag

- Multi compartment fold out bum bag
- Velcro securing to allow extra space for supplies
- Durable, water proof, light weight material
- Abrasion resistant material
- Three zip pouch window compartments
- Two zip side compartments
- Wide adjustable buckle belt clip
- Reflective strip



PR4	Fold Out Bum Bag - Red 270 x 160 x 45mm	A
-----	---	---

PR5	Fold Out Bum Bag - Green 270 x 160 x 45mm	A
-----	---	---

Training Manikins



Professional CPR/AED Training Manikins

Realistic to the eye and the touch and available individually or in convenient multipacks of four for class training purposes; these manikins are uniquely designed as a clamshell that accommodates an easy-to-insert face shield lung bag.

- The CPR Rate Monitor allows for instant feedback to both instructor and student regarding the rate of chest compression. The face/head tilt feature is unique and the manikin simulates the way an actual victim's head would move if he required CPR
- Realistic looking and acting manikin
- Easy to transport and set up
- Easy to insert face shield and lung bag
- Face shield is secured, no need for adhesives
- Full chest rise and lung bag ensures students receive the most accurate and realistic training by observing a true chest rise
- Face/head tilt demonstrates how to open the airway during CPR
- Clicker mechanism demonstrates how to administer adequate compressions

99700	Prestan Professional CPR/AED Manikin Infant	D
99700LB	Prestan Professional CPR/AED Lung Bags Infant	D
99700FS	Prestan Professional CPR/AED Face Shields Infant	D
99702	Prestan Professional CPR/AED Manikin Child	D
99702LB	Prestan Professional CPR/AED Lung Bags Child	D
99702FS	Prestan Professional CPR/AED Face Shields Child	D
99704	Prestan Professional CPR/AED Manikin Adult	D
99704LB	Prestan Professional CPR/AED Lung Bags Adult	D
99704FS	Prestan Professional CPR/AED Face Shields Adult	D
99706	Prestan Professional CPR/AED Manikin Family 2xAd1xCh1xInf	D

Cotec Finger Protection Tape



Cotec Safety Tape

Protects fingers, hands and wrists from minor injuries and abrasions in industrial processes.

Cotec tape can be used in any industrial operation where unprotected skin can be damaged. It can be applied selectively so that only the vulnerable area needs to be covered, reducing the cost while improving operator comfort.

The tape sticks to itself and does not leave any residue.

- Prevents cuts and sore fingers
- Protection against hot parts or mouldings
- Prevents finger prints on glass or polished surfaces
- Prevents product contamination by body oils and acids on fingers
- Guards against thread burns
- Quick to apply
- Open gauze allows skin to breathe
- Applying further layers increases protection

Cotec Safety Tape is available in 300mm packs cut into 25mm widths. Each roll is over 27 metres long.



In Industry

Assembly: To help protect fingers while pig-tailing, twisting and winding wires, or handling sharp components, fittings and circuit boards.

Hand sewing: For comfortable protection of fingers from abrasions during the sewing process.

Buffing and Grinding: Protect from abrasions in many machining and welding operations.

Mouldings: Protects the fingers when handling hot parts in moulding and pressing operations when gloves are not necessary.

80136	Cotec Finger Protection Tape Green - Pack 12	A
80137	Cotec Finger Protection Tape Blue - Pack 12	A

Disposable Gloves

MediKit's range of disposable gloves are:

- Powder Free: Low particle count reduces risk of product contamination
- Grip: Textured fingertips for improved grip
- Rolled Cuff: For additional strength and prevention of liquid roll back
- Freedom from Holes: AQL 1.5


Finex Latex Disposable Gloves

- Strong: Made from natural latex
- Sensitive Skin: Special processing reduces proteins and toxins which can cause sensitisation and irritation



40923	Finex Small - Box 100	B
40924	Finex Medium - Box 100	B
40925	Finex Large - Box 100	B
40926	Finex X-Large - Box 100	B

Sterile Gloves



40940	Sterile Latex Gloves - Pair	B
40941	Sterile Nitrile Gloves - Pair	B

Finite Powder Free Nitrile Disposable Gloves


- Strength: Synthetic rubber offers better puncture and abrasion resistance than natural rubber or vinyl
- Comfort: Soft elastic material feels cooler and is more comfortable and flexible to wear
- Sensitive Skin: No latex proteins – eliminates protein sensitisation




40942	Finite PF Medium - Box 100	B
40944	Finite PF Large - Box 100	B
40945	Finite PF X-Large - Box 100	B
40946	Finite PF XX-Large - Box 100	B

Finesse Powder Free Vinyl Disposable Gloves

- Sensitive: A thin powder free vinyl glove, with an excellent fit that feels like a second skin
- Sensitive Skin: No latex proteins - eliminates protein sensitisation
- Comfort: The soft vinyl formulation of Finesse PF provides greater comfort and is more flexible and dexterous to wear



40927	Finesse PF Medium - Box 100	B
40928	Finesse PF Large - Box 100	B
40929	Finesse PF X-Large - Box 100	B



Signs

SA711T	Pre-School First Aid Poster - A2	B
SA604T	Resuscitation of Children Poster - A2	B
SA598T	First Aid for Children Laminated Poster - A2	B
SA602T	First Aid for Eyes Laminated Poster - A2	B
SA603T	First Aid for Burns Laminated Poster - A2	B
SA597T	Safe Manual Handling Guidance Poster - A2	B
SA704T	COSHH Laminated Poster - A2	B
SA601T	Electric Shock Treatment Guide Poster - A2	B
SA660T	Anaphylaxis In Schools Poster - A2	B
SA3016	HSE Health & Safety Law Poster - A2	B
SA724T	Asthma Advice Poster - A2	B
SA9001	Snap Frame 25mm Silver - A2	B

